



Determination of veterinary drugs in foods of animal origin by QuEChERS coupled with ultra performance liquid chromatography-tandem mass spectrometry (UPLC-MS/MS)

Qianran Sun^a, Jun Liu^{a,*}, Yuan Gou^a, Tiejuan Chen^a, Xiaofang Shen^b, Tao Wang^a, Yongli Li^a, Huizhen He^a, Huidan Deng^a, Yi Hua^a

^a Technology Center, Chengdu Customs, Chengdu 610041, China

^b School of Food Science and Technology, Jiangnan University, Wuxi, 214122, China

ARTICLE INFO

Keywords:

Veterinary drugs
Foods of animal origin
UPLC-MS/MS
QuEChERS

ABSTRACT

A method using QuEChERS coupled with ultra-performance liquid chromatography-tandem mass spectrometry (UPLC-MS/MS) was developed for the determination of the residues of 19 veterinary drugs in ten animal-derived matrices, including beef, pork, sheep, horse, chicken, prawn, fish, liver, milk, and fat. This method was based on the enactment of veterinary drug compounds by Korea, Canada, the United States, and the European Union in recent years. The samples were extracted using 85% acetonitrile and separated on an ACQUITY UPLC HSS T3 column (2.1 mm × 100 mm, 1.8 μm) with a gradient elution of methanol-0.2% formic acid water as the mobile phase. The detection of the analytes was achieved through the use of positive ion electrospray ionization (ESI) and multiple reaction monitoring (MRM) modes, while the quantification was conducted via the matrix-matched external standard method. Following optimization, the linearity of the target veterinary residues in the ten matrices was observed to be satisfactory, having a range of 0.5–50.0 ng/mL ($R^2 > 0.991$). The limits of detection (LOD) were in the range of 0.01–1.29 μg/kg, while the limits of quantification (LOQ) were in the range of 0.02–4.31 μg/kg. The recoveries were observed to be in the range of 60.6–117.7 %, with relative standard deviations (RSDs) of ≤20.6 %. The method is straightforward and highly sensitive, and it satisfies the maximum limits set by the relevant standards of Korea, Canada, the USA, and the EU. It is well-suited for the rapid screening, qualitative, and quantitative analyses of metomidate, acetanilide, dl-methylephedrine, and other substances in foods of animal origin, providing technical assistance for cross-border food safety and testing.

1. Introduction

As the global demand for foodstuffs derived from animals continues to escalate, the concern surrounding the residues of veterinary drugs also rises. The presence of drug residues in food is typically attributed to the use of veterinary drugs in animal husbandry. In recent times, there have been a number of instances where veterinary drugs have been misused and abused, with the intention of enhancing the growth and feed efficiency of livestock [1]. Veterinary drug residues are defined as drug prototypes, metabolites, and drug impurities of toxicological significance. These residues accumulate or are stored in the edible parts of the food obtained from animals, including cells, tissues, or organs, or that enter the eggs of laying poultry or the milk of lactating animals, following the administration of veterinary drugs (including drug

additives). Prolonged ingestion can result in toxic effects on the human body [2], including carcinogenesis, organ dysfunction, drug resistance, and other adverse consequences [3–5]. To prevent these outcomes, pertinent regulations and standards regarding the maximum residue limit (MRL) of veterinary drugs have been promulgated in numerous countries (European, [6]).

The ongoing development and approval of novel veterinary pharmaceuticals has led to an expansion in the range of veterinary drugs present in foodstuffs. Consequently, the MRL for these drugs are subject to frequent updates and dynamic changes. To illustrate, drugs such as acetanilide, pentetrazol, methyl ephedrine hydrochloride, antipyrine, and guaifenesin are listed in the Korean Positive List with the corresponding specified limits in the respective animal-derived matrices. However, this is not the case in China, the European Union (EU),

* Corresponding author.

E-mail address: liujuncdhg@163.com (J. Liu).

<https://doi.org/10.1016/j.chroma.2025.465726>









Received 23 October 2024; Received in revised form 23 January 2025; Accepted 24 January 2025

Available online 26 January 2025

0021-9673/© 2025 Elsevier B.V. All rights are reserved, including those for text and data mining, AI training, and similar technologies.











Table 1

CAS number, molecular weight, solvent, chemical structure, and pharmacologically use of 19 veterinary residues.

NO.	Compound	CAS	Formula	Purity(%)	Solvent	Molecular weight	Chemical structure	Pharmacologically use
1	Acetanilide	103-84-4	C ₈ H ₉ NO	99.9	Acetonitrile	135.06		Aniline derivatives with analgesic and antipyretic properties
2	o-Aminobenzoic acid (Anthranilic acid)	118-92-3	C ₇ H ₇ NO ₂	99.9	Methanol	137.04		Intermediate for the synthesis of various drugs such as anti-inflammatory drugs and hypnotic drugs
3	Pentetrazol (Pentylene-tetrazol)	54-95-5	C ₆ H ₁₀ N ₄	99.9	Acetonitrile	138.09		A circulatory and respiratory stimulant
4	Methyl ephedrine hydrochloride (dl-Methylephedrine HCl)	18760-80-0	C ₁₁ H ₁₇ NO* HCl	99.9	Methanol	215.72		A non-selective adrenergic receptor agonist
5	Phenacetin	62-44-2	C ₁₀ H ₁₃ NO ₂	99.9	Methanol	179.09		An analgesic and fever-reducing drug
6	Antipyrine (Phenazone)	60-80-0	C ₁₁ H ₁₂ N ₂ O	98.0	Methanol	188.09		A nonsteroidal anti-inflammatory drug
7	Guaifenesin	93-14-1	C ₁₀ H ₁₄ O ₄	99.3	Methanol	198.02		Expectorants
8	Diethylcarbamazine	90-89-1	C ₁₀ H ₂₁ N ₃ O	99.8	Methanol	199.16		Treatment of certain worm infections



(continued on next page)

Table 1 (continued)

9	*Acriflavine(8048-52-0)	3,6-diamino-10-methylacridinium chloride	86-40-8	$C_{14}H_{14}N_3Cl$	99.9	Methanol	259.73		An antiseptic agent
		3,6-Acridinediamine(Proflavine)	92-62-6	$C_{13}H_{11}N_3$		Methanol	209.25		
10	Metomidate		5377-20-8	$C_{13}H_{14}N_2O_2$	99.0	Methanol	230.1		A sedative-hypnotic drug
11	Ethacridine lactate monohydrate (Acrinol)		6402-23-9	$C_{15}H_{15}N_3O^*$ $C_3H_6O_3^*H_2O$	94.6	Methanol	361.39		An antiseptic agent
12	Tripeleminamine hydrochloride		154-69-8	$C_{16}H_{21}N_3^*HCl$	99.9	Methanol	291.82		An antipruritic and first-generation antihistamine
13	Ormetoprim		6981-18-6	$C_{14}H_{18}N_4O_2$	99.9	Methanol	274.14		An antimicrobial agent
14	Sulfaethoxypyridazine		963-14-4	$C_{12}H_{14}N_4O_3S$	99.0	Acetonitrile	294.07		An antimicrobial agent
15	Yohimbine hydrochloride		65-19-0	$C_{21}H_{26}N_2O_3^*HCl$	97.7	Methanol	390.91		An alpha 2 adrenergic receptor antagonist
16	Buquinolate		5486-03-3	$C_{20}H_{27}NO_5$	98.3	Methanol	361.43		An anticoccidial agent
17	Dehydrocholic acid		81-23-2	$C_{24}H_{34}O_5$	99.6	Acetonitrile	402.24		A synthetic bile acid, hydrocholeretic

(continued on next page)

Table 1 (continued)

18	Loperamide hydrochloride	34552-83-5	$C_{29}H_{33}ClN_2 \cdot O_2 \cdot HCl$	98.7	Methanol	513.5		An antidiarrheal agent
19	Ciclesonide	126544-47-6	$C_{32}H_{44}O_7$	99.3	Acetonitrile	540.3		A glucocorticoid used to treat asthma and allergic rhinitis

*Commercially purchased acriflavine is a mixture of 3,6-diamino-10-methylacridinium chloride and 3,6-Acridinediamine (Proflavine).

Canada. Limits have been established for tripeleennamine hydrochloride in foods of animal origin in Korea and the United States (US). However, these limits are not mentioned in the Chinese standards GB 31,650–2019, GB 31,650.1–2022, and Announcement No. 250 of the Ministry of Agriculture and Rural Affairs. Furthermore, there is a dearth of corresponding detection methods for these drug residues. The abovementioned limits are also not covered in Canada and the EU. The limits for ormetoprim and buquinolate are included in Canada's

Maximum Residue Limits (MRLs) for veterinary drugs in foods. However, these limits have not been established in Korea, the EU, or the US. Ciclesonide is only referenced in the EU Regulation (EU) No 43/2020, and it has not been identified in other countries. In order to enhance the efficacy of veterinary drug testing standards and facilitate the acceleration of global trade, it is imperative to address the technical barriers to trade encountered by enterprises engaged in export activities. In this context, the standards of the corresponding veterinary drug detection

Table 2

Retention times and optimum operating mass spectrometric parameters for analysis of target veterinary drugs.

NO.	Analyte	Retention time	Ion source	Adduct ion composition	Precursor ion (m/z)	Product ion (m/z)	Cone voltage (V)	Collision energy (eV)
1	Diethylcarbamazine	4.94	ESI+	$[M + H]^+$	200.07	100.07*	2	14
					200.07	72.06	2	22
2	Methyl ephedrine hydrochloride	5.61	ESI+	$[M-HCl+H]^+$	180.10	147.04*	4	16
					180.10	117.07	4	18
3	Pentetrazol	6.17	ESI+	$[M + H]^+$	138.99	69.08*	28	22
					138.99	55.10	28	20
4	Ormetoprim	6.47	ESI+	$[M + H]^+$	275.03	123.08*	58	24
					275.03	81.03	58	42
5	Tripeleennamine hydrochloride	6.99	ESI+	$[M-HCl+H]^+$	256.05	119.18*	2	32
					256.05	91.04	2	30
6	Antipyrine	7.07	ESI+	$[M + H]^+$	189.01	56.08*	2	22
					189.01	104.06	2	26
7	o-Aminobenzoic acid	7.13	ESI+	$[M + H]^+$	137.94	92.07*	12	20
					137.94	65.05	8	24
8	3,6-Acridinediamine	7.14	ESI+	$[M + H]^+$	209.96	193.03*	44	30
					209.96	182.02	44	30
9	Yohimbine hydrochloride	7.25	ESI+	$[M-HCl+H]^+$	355.03	144.05*	6	28
					355.03	212.11	6	22
10	Acetanilide	7.37	ESI+	$[M + H]^+$	135.96	77.02*	18	20
					135.96	94.06	18	22
11	Sulfaethoxypyridazine	7.45	ESI+	$[M + H]^+$	295.08	156.00*	30	19
					295.08	140.00	30	20
12	Guaifenesin	7.52	ESI+	$[M + H]^+$	198.99	125.03*	28	8
					198.99	163.06	28	6
13	3,6-diamino-10-methylacridinium chloride	7.65	ESI+	$[M-HCl+H]^+$	223.91	182.02*	20	36
					223.91	209.05	20	30
14	Phenacetin	8.28	ESI+	$[M + H]^+$	180.10	110.20*	30	20
					180.10	138.00	30	30
15	Ethacridine lactate monohydrate	8.45	ESI+	$[M-C_3H_6O_3-H_2O+H]^+$	253.96	197.11*	20	30
					253.96	225.97	20	34
16	Metomidate	8.56	ESI+	$[M + H]^+$	230.97	126.99*	2	20
					230.97	95.01	2	8
17	Loperamide hydrochloride	9.31	ESI+	$[M-HCl+H]^+$	477.23	266.20*	50	25
					477.23	210.00	50	45
18	Dehydrocholic acid	9.67	ESI+	$[M + H]^+$	403.08	349.22*	18	20
					403.08	367.23	18	16
19	Buquinolate	11.3	ESI+	$[M + H]^+$	362.02	204.00*	10	38
					362.02	260.04	10	28
20	Ciclesonide	12.9	ESI+	$[M + H]^+$	541.31	323.20*	30	10
					541.31	523.20	30	20

* Monitored transition ions for quantification.




Fig. 1. The 19 target veterinary drug residues were separated effectively on the chromatogram (50ng/ml).

a. Diethylcarbamazine; b. Methyl ephedrine hydrochloride; c. Pentetrazol; d. Ormetoprim; e. Tripeleminamine hydrochloride; f. Antipyrine; g. o-Aminobenzoic acid; h-1. 3,6-Acridinediamine; i. Yohimbine hydrochloride; j. Acetanilide; k. Sulfaethoxypyridazine; l. Guaifenesin; h-2. 3,6-diamino-10-methylacridinium chloride; m. Phenacetin; n. Ethacridine lactate monohydrate; o. Metomidate; p. Loperamide hydrochloride; q. Dehydrocholic acid; r. Buquinolate; s. Ciclesonide

methods need to be instantly updated and improved simultaneously. Based on the differences in the regulatory limits between different countries, it is imperative to establish a universal analysis method for the residues of these veterinary drugs which have established maximum residue limits in only some countries and lack matching detection method standards. In this regard, Chae et al. [7] conducted a study on the detection of dl-methylephedrine hydrochloride in porcine muscle using LC-MS/MS. A study on the detection of flumethasone, dl-methylephedrine, and 2-hydroxy-4,6-dimethylpyrimidine in porcine muscle and pasteurised cow milk using the LC-MS/MS method was conducted by Zhang et al. [8]. Compared to this study, the above-mentioned investigations detected fewer veterinary substances, and the justification of the matrix was not as comprehensive. This study is therefore prospective in terms of both the types of veterinary drugs and the matrices tested. To the best of the authors' knowledge, no studies have been published currently on the use of high-performance liquid chromatography-tandem mass spectrometry (HPLC-MS/MS) methods for the detection of ciclesonide residues in foods of animal origin. Further comparative analyses with other relevant literature were presented in the comparison section(3.8) below.

A majority of pre-treatment purification methods for veterinary residues are based on solid-phase extraction (SPE), dispersive solid-phase extraction (DSPE), magnetic-solid phase extraction (M-SPE), dispersive liquid-liquid microextraction (DLLME), and liquid-liquid extraction (LLE) techniques [9–11]. At present, the utilisation of solid-phase extraction columns is relatively widespread. However, the series of operations, such as column activation, loading, washing, and elution, used in this pretreatment method are time-consuming and require a large number of reagents. Therefore, QuEChERS-DSPE was proposed as a pretreatment clean-up material [10]. This method requires straightforward operation, involves less solvent, effectively prevents sieve plate clogging due to complex matrices, and enhances the efficiency of pretreatment and purification of large-volume samples. The study involved ten matrices, namely beef, chicken, pork, prawn, fish, etc. The influencing factors, such as chromatographic conditions,

column selection, extraction solvent, purification powder, and centrifugation temperature, were optimized. The evaluation of specificity, linearity, matrix effect, limit of detection (LOD), limit of quantification (LOQ), recovery, and precision was employed as the basis for establishing a QuEChERS DSPE-UPLC-MS/MS method for the determination of 19 veterinary drug residues in foods of animal origin. This method provides a basis for the safety detection and market supervision of foods of animal origin and offers technical support for strengthening cross-boundary food safety assurance. The method is efficient, accurate, economical, and suitable for the rapid screening, confirmation, and quantitative detection of a wide range of matrices in animal-derived foods.

2. Material and methods

2.1. Reagents and chemicals

Methanol (MeOH) and acetonitrile (ACN) of LC-MS grade were purchased from CNW Technologies GmbH (Dusseldorf, Germany). Ammonium acetate ($\text{CH}_3\text{COONH}_4$), ammonium formate (HCOONH_4), formic acid (HCOOH), and ammonium hydroxide ($\text{NH}_3\cdot\text{H}_2\text{O}$) were of HPLC grade, obtained from Chengdu Kelong Chemical Co. Ltd. (Sichuan, China). Analytical-grade sodium chloride (NaCl), C18 adsorbent (40–63 μm , C18), anhydrous sodium sulfate (Na_2SO_4) and anhydrous magnesium sulphate (MgSO_4) were purchased from Anpel Laboratory Technologies (Shanghai, China). Additionally, the water used in the analyses was obtained from a Mili-Q Integral 3 ultrapure water device manufactured by Millipore Corporation (Boston, USA).

Buquinolate (10mg, 98.3%, CAS 5486–03–3), and Ciclesonide (100 mg, 99.3%, CAS 126,544–47–6) were procured from Beijing manhage bio-technology company (Beijing, China). All the other standards were bio-liquid reference materials with a purity of $\geq 94.6\%$, sourced from Alta scientific co. ltd. (Tianjin, China). Further details regarding the target veterinary residues, including the chemical abstract service (CAS) number, molecular weight, solvent, chemical structure, and




Fig. 2. Effect of optimising the choice of extraction solvents on the recovery of 19 veterinary drugs.




Fig. 2. (continued).

pharmacologically use are provided in Table 1.

Most notably, acriflavine is a known mixture of the parent compound (3,6-diamino-10-methylacridinium chloride) and the demethylated derivative, proflavine, as indicated by the purchased standard [12,13]. The concentration ratio of chlorinated 3,6-diamino-10-methylacridine to 3,6-acridinediamine was 2:1. The chromatographic columns exhibited distinct peaks and retention times. Consequently, these substances were treated as two separate substances in the following analyses.

2.2. Standard solutions

2.2.1. Preparation of standard stock solution

0.01 g of ciclesonide and buquinolate were weighed accurately in a 10 mL volumetric flask. Acetonitrile and methanol were added to dissolve and dilute the standard to the scale, respectively. The individual standard stock solution was prepared at a mass concentration of 1000 µg/mL and stored at -20°C , protected from light. The remaining liquid

standards were accurately pipetted at 1 mL each, with the addition of methanol or acetonitrile to achieve a final volume of 10 mL. The standard stock solutions were then prepared with a mass concentration of either 10 µg/mL or 100 µg/mL, and stored at -20°C in a dark environment for future use. An appropriate amount of the standard stock solutions of the abovementioned veterinary residues was accurately measured, diluted with methanol, and formulated into a mixed standard working stock solutions with a mass concentration of 1 µg/mL. Subsequently, a series of working standard solutions were formulated at concentrations of 0.5, 1.0, 2.0, 3.0, 5.0, 10.0, 20.0, and 50.0 ng/mL, for UPLC-MS/MS detection.

2.2.2. Preparation of blank matrix curves

The mixed standard working solution was quantified and incorporated into the residues of eight extracted and purified blank samples (the pretreatment method was identical to that described in Section 2.3), and then blown to dryness under nitrogen at 50°C . 1 mL of methanol was




Fig. 3. Recoveries of 19 veterinary residues at different optimized centrifugation temperatures (take grass carp for example).

added to the resulting samples, followed by vortexing and dissolution. The resulting solutions were then synthesized to form a matrix-matched series of mixed standard solutions at concentrations of 0.5, 1.0, 2.0, 3.0, 5.0, 10.0, 20.0, and 50.0 ng/mL. The solutions were poured through a 0.22- μ m filter membrane for UPLC-MS/MS on-line determination. The standard curve was plotted with the measured characteristic ion peak area as the vertical coordinate and the corresponding standard solution concentration as the horizontal coordinate.

2.3. Sample preparation

A total of 50 or more batches of muscle samples of cattle (Muscle), swine (Muscle), sheep (Muscle), horse (Muscle), chicken (Muscle), fish (Grass carp, Skin and Muscle), prawn (Muscle), liver (pork liver), milk (cow milk), and fat (lard oil) were used for the experiments. All the samples were obtained from local markets or supermarkets. Each sample was homogenized in a homogenizer and stored at -20°C .

2.3.1. Beef, sheep, horse, chicken, fish, prawn, liver, milk, and fat samples

2.0 g of the sample was accurately weighed and placed individually into 50-ml screw-capped centrifuge tubes. Subsequently, the tubes were spiked with an appropriate quantity of mixed standard solutions and kept in 4°C overnight. This process enables the solution to penetrate completely into the tissue cells (the fat samples were heated at 60°C in a water bath until they melted). A volume of 10 mL of 85% acetonitrile was added to the sample (the addition of 1 g NaCl for milk samples), which was then extracted by shaking for 5 mins. The sample was then extracted by an ultrasonic water bath for 15 mins, after which it was centrifuged at 10,000 r/min at 4°C for 5 mins.

The supernatant was added to a purification tube containing 100 mg of C18 and 1100 mg of anhydrous magnesium sulfate, vortexed and mixed for 2 mins, and then subjected to freeze-centrifugation at 10,000 r/min and 4°C for 5 mins. The final supernatant was transferred to a 15 mL glass nitrogen blowing tube and blown to dryness at 45°C . Finally, 1 mL of methanol was accurately added, vortexed to dissolve the residue, and filtered through a 0.22- μ m microporous membrane for UPLC-MS/MS on-line analysis.

2.3.2. Pork samples

2.0 g of the sample was accurately weighed and placed individually into 50-ml screw-capped centrifuge tubes. Subsequently, the tubes were

spiked with an appropriate quantity of mixed standard solutions and allowed to stand for 30 mins at room temperature, which enables the solution to penetrate completely into the tissue cells. A volume of 10 mL of 85% acetonitrile was added to the sample, which was then extracted by shaking for 5 mins. The sample was subsequently extracted by an ultrasonic water bath for 15 mins, after which it was centrifuged at 10,000 r/min at 4°C for 5 mins. The supernatant was combined with 5 mL of hexane saturated with acetonitrile, vortexed for 2 min, and then freeze-centrifuged at 10,000 r/min 4°C for 5 mins.

The acetonitrile layer was transferred to a purification tube containing 100 mg of C18 and 1100 mg of anhydrous magnesium sulfate, vortexed for 2 min, and then freeze-centrifuged at 10,000 r/min 4°C for 5 min. The final supernatant was pipetted into a 15-mL glass nitrogen blowing tube and nitrogen-blown to a clean dry state at 45°C . 1 mL of methanol was added accurately to dissolve the residue and then passed through a 0.22- μ m filter membrane for UPLC-MS/MS analysis.

2.4. Instrumentation

2.4.1. Chromatographic conditions

The detection of 19 veterinary drug residues was performed on a UPLC system (UPLC I-class-XEVOTQ-XS, Waters Corporation, USA). The chromatographic separation was performed on an Acquity UPLC HSS T3 (2.1 mm $\times$ 100 mm, 1.8 μ m, Waters Corporation, USA).

Methanol (mobile phase A) was used as the organic phase, and 0.2% formic acid (HCOOH) in distilled water (mobile phase B) was used as the aqueous phase. The optimal elution conditions for the mixture of veterinary drug residues used were 5 % A at 0–2 min, 5–90 % A at 2–10 min, 90–95 % A at 10–10.1 min, 95 % A at 10.1–12 min, 95–5 % A at 12–12.5 min, and 5 % A at 12.5–15 min. The total duration of the experiment was 15 mins. Moreover, the flow rate was 0.25 mL/min, and the injection volume was 3.0 μ L.

2.4.2. Detection conditions

Mass spectrometry was conducted in the electrospray positive ion (ESI+) mode. The electrospray capillary voltage was set at 3.00 kV, while the ESI+ cone voltage was fixed at 30 V. The dissolution gas flow rate was adjusted to 1000 L/h, while the collision gas flow rate was set at 0.15 mL/min at a dissolution temperature of 500°C . Additionally, a qualitative analysis of each compound was conducted in the multiple reaction monitoring mode (MRM).




Fig. 4. Recoveries of veterinary residues at different types of purified powder (take liver for example). Group1: 80 mg C18+1000 mg MgSO₄+80 mg PSA; Group2: 80 mg C18+1000 mg MgSO₄+40 mg GCB; Group3: 80 mg C18+1000 mg MgSO₄+80 mg PSA+40 mg GCB; Group4: 80 mg C18+1000 mg MgSO₄.

The UPLC-MS/MS system was controlled by using the MassLynx™ software. Data were processed through the TargetLynx™ software (both from Waters) to obtain standard curves, linear equations, and calculated results. Microsoft Excel was employed to process the data and draw graphs. The spectra were plotted by using SciDAVis and Origin.

3. Results and discussion

3.1. Optimisation of detection conditions

A single standard solution (100 ng/mL) was selected for the experiment to determine the precursor ions of the analyte in both positive and negative ion modes. Based on the optimised detector voltage conditions, the two highest-responding fragment ions were selected as their respective product ions. As per a previous study, the most common ion adducts observed in positive ion electrospray ionization (ESI+) mass spectrometry are $[M + H]^+$, $[M+NH_4]^+$, $[M+Na]^+$, $[M + K]^+$, $[2M+H]^+$, $[2M+NH_4]^+$, $[2M+Na]^+$ [14]. Ionised acidic residues (Sulfonylamino, carboxyl, phenolic hydroxyl, etc.) readily form ion pairs with positive ions such as Na^+ and NH_4^+ [15–17]. The results revealed that most of the target veterinary drug residues formed $[M + H]^+$ adducts (Table 2), with higher responses in the positive ion mode. The selection of the positive ion mode for the detection of substances with similar aniline or amino ($-NH_2$) structures was reported by Palm et al. [18], and Mezghich et al. [19]. The optimal mass spectrometry parameters of each target substance are presented in Table 2.

3.2. Chromatographic column selection

The separation of 19 target veterinary residues on ACQUITY UPLC BEH C18 (100 mm × 2.1 mm, 1.7 μm, Waters, USA), ACQUITY UPLC HSS T3 (2.1 mm × 100 mm, 1.8 μm, Waters, USA), CAPCELL PAK C18 MGIII (2.0 mm × 100 mm, 5 μm, Shiseido, Japan), and InertSustain C18 (2.1 mm × 150 mm, 3 μm, GL Science, Japan) columns was compared. Although the ACQUITY UPLC BEH C18 column is commonly employed


as a general-purpose ultra-high performance chromatographic column for a diverse array of analytes, the separation of the target veterinary residues was found to be suboptimal. The compounds demonstrated comparable performance on both the InertSustain C18 and MGIII columns. Specifically, dehydrocholic acid and 3,6-acridinediamine presented elevated baselines on the InertSustain C18 column. In contrast, the HSS T3 column demonstrated a superior peak shape and response. The HSS T3 stationary phase is a C18 stationary phase that is compatible with a 100% aqueous mobile phase [20,21], which is more suitable for the retention of water-soluble and polar small organic compounds. A majority of the target veterinary residues possess a high degree of polarity. Walter et al. [22] employed a combination of six polar compounds to assess the retention and selectivity of distinct column stationary phases. The authors state that for the positively charged polar analytes, the HSS T3 column exhibited the most retention. Moreover, Medellín-Martínez et al. [23] used ACQUITY UPLC HSS T3 to separate clenbuterol residues from beef. Therefore, the T3 column was chosen as the analytical column for the experiment.

3.3. Optimisation of mobile phases

In general, polar solvents (e.g., water, methanol, acetonitrile, etc.) are employed in ESI-MS experiments owing to their propensity to undergo electrochemical reactions in the spraying nozzle [24]. The peak shape behaviors obtained using methanol and acetonitrile were compared. It was observed that the methanol/water mobile phase exhibited a lower elution capacity than the acetonitrile/water mobile phase. Consequently, the retention time of the target veterinary residues using the methanol/water mobile phase was prolonged, thereby facilitating the separation of the chromatographic peaks. Following a comprehensive evaluation, methanol was selected as the organic phase.

The separation effects and response discrepancy of different water phases for each target were investigated using 0.1% HCOOH, 0.2% HCOOH, 5 mmol/L HCOONH₄, 5 mmol/L CH₃COONH₄, and 5 mmol/L CH₃COONH₄ + 0.1% HCOOH as the mobile aqueous phase,

a



b




Fig. 5. Optimisation of C18 dosage for 19 veterinary residues in ten matrices.

respectively. The results demonstrated that the addition of HCOOH could effectively improve the peak shape and protonation ability, thereby facilitating the formation of the $[M + H]^+$ peak and augmenting the response [25]. It was observed that compared with 0.1% HCOOH, 0.2% HCOOH exhibited a higher response and good peak shape. Therefore, methanol-0.2% formic acid water was finally selected as the mobile phase. Under optimal conditions, the rapid and complete

separation of the 19 target veterinary drug residues was achieved on the chromatograms as shown in Fig. 1.

3.4. Optimization of sample extraction

3.4.1. Extraction solvent

Extraction reagents demonstrate a significant effect on the recovery,




Fig. 6. Optimisation of anhydrous magnesium sulfate (MgSO_4) for 19 veterinary residues in ten matrices.

the time of appearance of the target peaks, and the peak shape of several veterinary drug residues. Commonly used extraction solvents in foods of animal origin include methanol, acetonitrile, and ethyl acetate [26,27]. Furthermore, the acidity or alkalinity of the extraction environment influences the efficiency of the extraction process with regard to the target substances. In order to compare the extraction efficiency of the substances in question, the study employed formic acid and ammonium acetate as buffer solutions [28]. Acetonitrile, ethyl acetate, 0.1% formic acid acetonitrile, and 1% ammonia acetonitrile were chosen as extraction solvents to compare the extraction effect. The results presented in Fig. 2a reveal that acetonitrile serves as an effective extraction solvent in fish matrices, exhibiting notable benefits. The recoveries of metomidate and acrinol in acetonitrile and 0.1% formic acid acetonitrile were comparable, whereas sulfaethoxypyridazine exhibited a higher recovery in the latter solvent. In the case of prawn (Fig. 2b), o-aminobenzoic acid, 3,6-acridinediamine, and tripeleannamine were better recovered when 0.1% formic acid acetonitrile was used as the extraction solvent. However, guaifenesin, metomidate, buquinolate, ormetoprim, methyl ephedrine hydrochloride, and yohimbine were not suitable for extraction in acidic environments, with recoveries below 30%. The extractions of pentetrazol, diethylcarbamazine, dehydrocholic acid, and loperamide were more effective when ammonia acetonitrile was used, while the remaining substances demonstrated higher recoveries under acetonitrile conditions. In chicken (Fig. 2c), o-aminobenzoic acid, ormetoprim, and buquinolate showed an extremely high recovery efficiency (up to 120% or more) when extracted by 1% ammonia acetonitrile. On the other hand, acetanilide, methyl ephedrine hydrochloride, metomidate, antipyrine, and loperamide demonstrated lower recovery efficiencies (<50%). The recovery of all substances by using ethyl acetate was limited to a maximum value of 40 % by ethyl acetate. A majority of substances demonstrated a good recovery efficiency in acetonitrile. Similarly, acetonitrile demonstrated high extraction efficiencies in pork and beef matrices. However, such high values were not observed for buquinolate and dehydrocholic acid in pork (Fig. 2d), and diethylcarbamazine, 3,6-diamino-10-methylacridinium chloride, acrinol, and

buquinolate in beef matrices (Fig. 2e).

Ethyl acetate also exhibited the worst recovery efficiency, which was hardly higher than 20%. Given that animal matrices contain a large amount of lipids, and many lipophilic compounds are readily extracted by ethyl acetate, the degree of matrix interference is high, resulting in a low extraction efficiency [29]. The extreme emulsification of ethyl acetate as an extraction solvent is also another possible factor. Acetonitrile exhibits good solubility and strong penetration, which can cause precipitation of proteins and effectively avoid over-extraction of lipids [30, 31]. Therefore, it has the advantage of the ability to extract the target more effectively than other solvents in most matrices. However, proteins are subject to agglutination in organic solvents, and the recoveries of certain substances under acetonitrile extraction conditions were found to be unsatisfactory. The addition of an appropriate amount of purified water to the extraction reagent can reduce the proportion of acetonitrile and promotes sample dispersion, as well as decrease the rate of homogeneous agglomeration [32]. Furthermore, a mixture of acetonitrile and water is capable of extracting a wide polarity range of analytes from the matrix [33]. Therefore, the experiment was continued to compare the extraction efficiencies of 80%, 85%, 90% and 95% acetonitrile in Fig. 2. The results revealed that 85% acetonitrile exhibited good recovery efficiencies for the target substances in all the studied matrices, especially for buquinolate and dehydrocholic acid, with a significant enhancement effect. The subsequent optimized treatment could also improve the sensitivity of the method. Thus, 85% acetonitrile in water was finally chosen as the extraction solvent.

The high lipid content of pork necessitated the implementation of an additional extraction step. Animal fats are composed primarily of long-chain saturated fatty acids, which are distinguished by the presence of two groups: a polar carboxyl group and a non-polar hydrocarbon group. The solubility of the fatty acids in water is influenced by the length of the hydrocarbon group, with longer hydrocarbon groups exhibiting reduced polarity and increased non-polarity. The target veterinary residues are predominantly polar substances with a high polarity. Therefore, it is advisable to avoid over-extraction with highly polar solvents in order to

Table 3

Matrix effects of targeted analytes in ten matrices.

No.	Analyte	Solvent standard curve	Correlation coefficient(r^2)	Matrix	Matrix-matched calibration curve	Correlation coefficient(R^2)	Matrix effect (ME/%)
1	Acetanilide	$Y = 61,204.3X + 12,316.2$	0.998	Beef	$Y = 12,924.03X - 6334.18$	0.999	-78.88
				Pork	$Y = 29,574.1X - 7721.23$	0.992	-51.68
				Sheep	$Y = 55,623.2X + 16,368.8$	0.999	-9.12
				Horse	$Y = 52,837.3X + 22,144.1$	0.997	-13.67
				Chicken	$Y = 24,353.8X + 75,795.2$	0.995	-60.21
				Prawn	$Y = 7768.35X - 2955.34$	0.993	-87.31
				Fish	$Y = 19,834.3X + 1309.34$	0.998	-67.59
				Liver	$Y = 31,928.4X + 14,333.1$	0.994	-47.83
				Milk	$Y = 58,934.3X + 43,451.2$	0.999	-3.71
				Fat	$Y = 68,368.8X + 50,030.5$	0.999	11.71
2	o-Aminobenzoic acid	$Y = 12,189.2X + 8143$	0.996	Beef	$Y = 6983.77X + 73,455.6$	0.993	-42.71
				Pork	$Y = 7938.26X + 25,469.2$	0.993	-34.87
				Sheep	$Y = 11,356.3X + 234,554$	0.999	-6.83
				Horse	$Y = 11,981.1X + 96,009.3$	0.996	-1.71
				Chicken	$Y = 6425.02X + 18,739.3$	0.993	-47.29
				Prawn	$Y = 3483.35X + 9345.87$	0.994	-71.42
				Fish	$Y = 7309.32X + 39,285.6$	0.992	-40.03
				Liver	/	/	/
				Milk	$Y = 11,021.5X + 198,743$	0.995	-9.58
				Fat	$Y = 13,364.9X + 19,809$	0.998	9.65
3	Pentetrazol	$Y = 21,827.3X + 923.434$	0.997	Beef	$Y = 8234.43X + 8782.33$	0.999	-62.27
				Pork	$Y = 6543.33X - 532.384$	0.997	-70.02
				Sheep	$Y = 16,748.8X + 9833.29$	0.999	-23.27
				Horse	$Y = 14,537.4X + 18,384.3$	0.994	-33.40
				Chicken	$Y = 7473.45X + 38,921.3$	0.994	-65.76
				Prawn	$Y = 3974.87X + 12,242.5$	0.997	-81.79
				Fish	$Y = 5476.86X + 6393.44$	0.995	-74.91
				Liver	$Y = 7783.22X + 4894.36$	0.999	-64.34
				Milk	$Y = 18,739.2X - 3338.35$	0.999	-14.15
				Fat	$Y = 23,293.7X + 2344.54$	0.999	12.82
4	Methyl ephedrine hydrochloride	$Y = 26,538.4X + 2534.4$	0.997	Beef	$Y = 5283.59X + 9123.67$	0.998	-80.09
				Pork	$Y = 3928.5X + 2896.88$	0.996	-85.20
				Sheep	$Y = 18,274.8X - 2122.61$	0.999	-31.14
				Horse	$Y = 8763.08X + 672.34$	0.997	-66.98
				Chicken	$Y = 6473.23X - 1673.57$	0.996	-75.61
				Prawn	$Y = 7364.23X + 11,245.2$	0.993	-72.25
				Fish	$Y = 7683.7X - 2456.78$	0.999	-71.05
				Liver	$Y = 3789.71X - 1778.63$	0.993	-85.72
				Milk	$Y = 18,739.3X - 8123.89$	0.999	-29.39
				Fat	$Y = 39,893.5X - 2332.74$	0.999	50.32
5	Phenacetin	$Y = 108393X + 2631.71$	0.997	Beef	$Y = 12,783.1X - 14,238.1$	0.998	-88.21
				Pork	$Y = 17,849.4X - 8877.55$	0.997	-83.53
				Sheep	$Y = 89,837.1X - 2337.6$	0.999	-17.12
				Horse	$Y = 79,803.2X - 9762.95$	0.995	-26.38
				Chicken	$Y = 14,757.32X - 18,780.3$	0.995	-86.39
				Prawn	$Y = 9793.84X - 387.554$	0.994	-90.96
				Fish	$Y = 15,243.42X - 28,767.6$	0.992	-85.94
				Liver	$Y = 49830X + 10,983.1$	0.997	-54.03
				Milk	$Y = 68,973.7X - 7674.56$	0.999	-36.37
				Fat	$Y = 118238X - 4333.14$	0.999	9.08
6	Antipyrine	$Y = 48,372.1X + 1359.72$	0.998	Beef	$Y = 7832.43X + 1976.64$	0.997	-83.81
				Chicken	$Y = 33,243.9X + 56,752.9$	0.995	-31.27
				Pork	$Y = 46,453.3X - 8909.54$	0.994	-3.97
				Sheep	$Y = 44,354.1X + 7808.32$	0.999	-8.31
				Horse	$Y = 31,232.5X + 8797.73$	0.998	-35.43
				Prawn	$Y = 7342.65X + 2309.43$	0.997	-84.82
				Fish	$Y = 13,245.7X - 7909.6$	0.997	-72.62
				Liver	$Y = 27,242.2X + 5098.23$	0.993	-43.68
				Milk	$Y = 32,343.3X + 2654.89$	0.998	-33.14
				Fat	$Y = 54,675.5X - 4984.45$	0.998	13.03
7	Guaifenesin	$Y = 8435.31X + 2423.65$	0.998	Beef	$Y = 5865.13X - 209.453$	0.993	-30.47
				Pork	$Y = 1453.32X - 988.032$	0.996	-82.77
				Sheep	$Y = 7435.22X + 2988.65$	0.994	-11.86
				Horse	$Y = 6675.71X + 3767.77$	0.999	-20.86
				Chicken	$Y = 2149.42X + 9943.3$	0.999	-74.52
				Prawn	$Y = 7345.96X - 715.209$	0.999	-12.91
				Fish	$Y = 1764.43X - 2123.69$	0.994	-79.08
				Liver	$Y = 4324.6X - 885.232$	0.999	-48.73
				Milk	$Y = 7453.32X - 402.143$	0.999	-11.64
				Fat	$Y = 11,234.8X - 2541.06$	0.999	33.19
8	Diethylcarbamazine	$Y = 167543X + 32,468.9$	0.998	Beef	$Y = 25,432.8X + 178,633$	0.999	-84.82
				Pork	$Y = 37,643.2X + 3789.45$	0.995	-77.53
				Sheep	$Y = 45,643.4X - 5981.2$	0.999	-72.76

(continued on next page)

Table 3 (continued)

No.	Analyte	Solvent standard curve	Correlation coefficient(r^2)	Matrix	Matrix-matched calibration curve	Correlation coefficient(R^2)	Matrix effect (ME/%)
9	3,6-diamino-10-methylacridinium chloride	$Y = 21,928.9X + 134,761$	0.998	Horse	$Y = 59,869.9X + 17,867.5$	0.998	-64.27
				Chicken	$Y = 65,436.5X + 139,733$	0.999	-60.94
				Prawn	$Y = 75,355.7X + 208,945$	0.995	-55.02
				Fish	$Y = 37,549.1X + 8554.53$	0.993	-77.59
				Liver	$Y = 65,413.8X - 5234.61$	0.993	-60.96
				Milk	$Y = 134636X - 48,891.1$	0.999	-19.64
				Fat	$Y = 223,411.3X + 16,230.9$	0.998	33.35
				Beef	$Y = 4342.87X + 1323.21$	0.992	-80.20
				Pork	$Y = 7324.43X - 402.334$	0.992	-66.60
				Sheep	$Y = 19,876.7X - 2543.87$	0.999	-9.36
				Horse	$Y = 12,321.3X - 343.724$	0.997	-43.81
				Chicken	$Y = 7342.89X - 7342.72$	0.995	-66.52
				Prawn	$Y = 5342.23X + 7897.66$	0.998	-75.64
				Fish	$Y = 6342.95X + 3972.34$	0.993	-71.07
				Liver	$Y = 4843.99X + 5698.25$	0.995	-77.91
				Milk	$Y = 15,462.9X - 6988.53$	0.997	-29.49
10	3,6-Acridinediamine	$Y = 31,288.7X + 4341.08$	0.997	Fat	$Y = 23,343.2X - 4874.78$	0.998	6.45
				Beef	$Y = 5231.24X + 387.564$	0.999	-83.28
				Pork	$Y = 7651.32X + 3987.21$	0.998	-75.55
				Sheep	$Y = 24,533.6X - 1987.72$	0.999	-21.59
				Horse	$Y = 24,324.6X - 711.89$	0.994	-22.26
				Chicken	$Y = 15,343.2X + 38,973.1$	0.998	-50.96
				Prawn	$Y = 6451.16X + 4897.34$	0.999	-79.38
				Fish	$Y = 7923.67X - 5782.13$	0.998	-74.68
				Liver	$Y = 9874.46X + 3768.71$	0.998	-68.44
				Milk	$Y = 26,548.6X - 8972.19$	0.999	-15.15
				Fat	$Y = 36,433.2X + 9553.8$	0.999	16.44
				Beef	$Y = 23,255.2X - 49,807.9$	0.998	-87.94
				Pork	$Y = 48,984.3X - 46,578$	0.994	-74.60
				Sheep	$Y = 229,381.2X - 20,391.1$	0.999	18.95
				Horse	$Y = 218,307.2X - 11,135.2$	0.998	13.21
				Chicken	$Y = 38,324.1X - 19,875.3$	0.991	-80.13
11	Metomidate	$Y = 192,831.2X + 11,113.4$	0.998	Prawn	$Y = 24,342.2X + 28,795.9$	0.999	-87.38
				Fish	$Y = 21,732.7X - 34,098$	0.996	-88.73
				Liver	$Y = 132458X + 2223.23$	0.997	-31.31
				Milk	$Y = 167,434.3X - 57,686.3$	0.999	-13.17
				Fat	$Y = 255332X - 17,817.2$	0.998	32.41
				Beef	$Y = 26,732.2X + 737.092$	0.999	-83.12
				Pork	$Y = 40,342.6X - 6094.42$	0.996	-74.53
				Sheep	$Y = 130293X - 19,878.6$	0.999	-17.73
				Horse	$Y = 149303X - 2898.1$	0.998	-5.73
				Chicken	$Y = 33,442.6X - 26,664.3$	0.995	-78.88
				Prawn	$Y = 23,244.9X + 28,904$	0.993	-85.32
				Fish	$Y = 36,553.3X - 56,375.5$	0.994	-76.92
				Liver	$Y = 49,736.8X + 25,986.4$	0.996	-68.60
				Milk	$Y = 149280X - 48,277.8$	0.994	-5.74
				Fat	$Y = 198,374.1X - 17,866.6$	0.998	25.26
12	Ethacridine lactate monohydrate	$Y = 158373X + 23,879$	0.995	Beef	$Y = 47,638.4X + 3873.5$	0.997	-77.57
				Pork	$Y = 68,394.1X - 4034.3$	0.994	-67.79
				Sheep	$Y = 224,221.2X - 4164.3$	0.999	5.59
				Horse	$Y = 209238X - 10,989.8$	0.998	-1.47
				Chicken	$Y = 97,384.7X + 189,775$	0.997	-54.14
				Prawn	$Y = 41,211.8X + 49,811.2$	0.994	-80.59
				Fish	$Y = 58,927.5X - 34,532.7$	0.998	-72.25
				Liver	$Y = 49,383.7X + 28,934.3$	0.991	-76.74
				Milk	$Y = 183948X - 54,453$	0.999	-13.38
				Fat	$Y = 283,943.2X - 267.986$	0.996	33.71
				Beef	$Y = 10,096.23X + 21,133.7$	0.999	-51.76
				Pork	$Y = 14,533.3X - 1360.45$	0.998	-30.56
				Sheep	$Y = 8453.23X + 378.123$	0.997	-59.61
				Horse	$Y = 6782.43X + 354.100$	0.999	-67.59
				Chicken	$Y = 21,537.4X + 48,938.2$	0.998	2.91
14	Ormetoprim	$Y = 20,928.1X + 2132.16$	0.997	Prawn	$Y = 14,539.3X + 6549.7$	0.993	-30.53
				Fish	$Y = 16,758.8X + 4343.32$	0.999	-19.92
				Liver	$Y = 5340.64X + 694.122$	0.992	-74.48
				Milk	$Y = 20,246.1X - 5193.2$	0.996	-3.26
				Fat	$Y = 30,322.2X - 2123.46$	0.999	44.89
				Beef	$Y = 13,452.6X + 2389.37$	0.999	-72.41
				Pork	$Y = 17,654.4X - 7776.68$	0.998	-63.80
				Sheep	$Y = 53,420.4X - 8123.56$	0.997	9.55
				Horse	$Y = 50,424.2X - 3111.64$	0.998	3.40
				Chicken	$Y = 28,764.7X + 63,452$	0.992	-41.01
				Prawn	$Y = 9975.54X + 21,523.1$	0.996	-79.54
				Fish	$Y = 15,439.1X - 756.657$	0.999	-68.34
				Liver	$Y = 16,532.2X + 23,576.8$	0.998	-66.10
15	Sulfaethoxypyridazine	$Y = 48,764.1X + 1233.2$	0.992				

(continued on next page)

Table 3 (continued)

No.	Analyte	Solvent standard curve	Correlation coefficient(r^2)	Matrix	Matrix-matched calibration curve	Correlation coefficient(R^2)	Matrix effect (ME/%)
16	Yohimbine hydrochloride	$Y = 129301X + 6028.2$	0.992	Milk	$Y = 53,424.5X - 25,779.9$	0.999	9.56
				Fat	$Y = 79,342.5X + 27,774.5$	0.998	62.71
				Beef	$Y = 24,938.5X + 21,187.3$	0.997	-80.71
				Pork	$Y = 38,971.4X + 5435.2$	0.998	-69.86
				Sheep	$Y = 88,394.1X - 15,647.2$	0.998	-31.64
				Horse	$Y = 76,473.2X - 4988.52$	0.998	-40.86
				Chicken	$Y = 78,326.9X + 62,222$	0.999	-39.42
				Prawn	$Y = 21,122.4X + 112,435$	0.998	-83.66
				Fish	$Y = 33,242.4X - 1879.23$	0.993	-74.29
				Liver	$Y = 21,113.6X + 14,561.3$	0.993	-83.67
				Milk	$Y = 98,337.6X - 49,876.9$	0.995	-23.95
				Fat	$Y = 133459X + 2777.25$	0.999	3.22
17	Buquinolate	$Y = 117865X + 7686.45$	0.993	Beef	$Y = 13,242.4X - 189.335$	0.995	-88.76
				Pork	$Y = 63,243.6X + 18,971.6$	0.992	-46.34
				Sheep	$Y = 155432X - 33,457.1$	0.999	31.87
				Horse	$Y = 85,334.7X - 9878.41$	0.996	-27.60
				Chicken	$Y = 41,234.2X + 18,666.8$	0.997	-65.02
				Prawn	$Y = 34,345.1X + 15,779.3$	0.999	-70.86
				Fish	$Y = 64,242.1X - 1223.12$	0.997	-45.50
				Liver	$Y = 58,758.7X + 7534.88$	0.996	-50.15
				Milk	$Y = 104246X - 31,231.2$	0.999	-11.55
				Fat	$Y = 135334X + 564.735$	0.999	14.82
				Beef	$Y = 3346.14X - 278.712$	0.995	-18.85
				Pork	$Y = 1768.31X - 454.567$	0.996	-57.11
18	Dehydrocholic acid	$Y = 4123.21X + 802.332$	0.994	Sheep	$Y = 4340.11X - 945.762$	0.999	5.26
				Horse	$Y = 3976.34X - 1564.67$	0.998	-3.56
				Chicken	$Y = 2568.37X + 3432.08$	0.993	-37.71
				Prawn	$Y = 5785.37X + 567.098$	0.997	40.31
				Fish	$Y = 2345.15X - 218.225$	0.994	-43.12
				Liver	$Y = 2123.21X - 543.334$	0.999	-48.51
				Milk	$Y = 2531.11X - 521.567$	0.991	-38.61
				Fat	$Y = 4345.45X - 967.204$	0.997	5.39
				Beef	$Y = 50,962.1X - 2112.34$	0.999	-85.62
				Pork	$Y = 84,532.2X - 34,347.3$	0.993	-76.14
				Sheep	$Y = 414243X - 64,341.5$	0.999	16.90
				Horse	$Y = 371461X - 1435.7$	0.998	4.83
19	Loperamide hydrochloride	$Y = 354353X + 43,516.5$	0.997	Chicken	$Y = 155322X + 245,609$	0.998	-56.17
				Prawn	$Y = 55,740.1X + 187,675$	0.998	-84.27
				Fish	$Y = 49,234.1X - 69,821.9$	0.994	-86.11
				Liver	$Y = 193241X + 48,977.5$	0.997	-45.47
				Milk	$Y = 265309X - 91,123$	0.993	-25.13
				Fat	$Y = 425632X + 89,856.4$	0.999	20.12
				Beef	$Y = 4532.42X - 309.819$	0.993	-57.46
				Pork	$Y = 6342.17X - 235.345$	0.996	-40.47
				Sheep	$Y = 4344.23X - 563.074$	0.999	-59.23
				Horse	$Y = 4453.22X - 193.321$	0.997	-58.20
				Chicken	$Y = 7543.64X - 112.133$	0.992	-29.20
				Prawn	$Y = 4764.21X - 242.456$	0.997	-55.28
20	Ciclesonide	$Y = 10,654.3X - 6236.54$	0.994	Fish	$Y = 3985.23X - 135.411$	0.996	-62.60
				Liver	/	/	/
				Milk	$Y = 8982.1X - 10,788.2$	0.993	-15.70
				Fat	$Y = 6563.31X + 1358.46$	0.998	-38.40

minimize losses. According to the similarity compatibility principle, nonpolar organic solvents interact with the long hydrophobic chains of fatty acids and neutral lipids through hydrophobic interactions (van der Waals forces), thereby solubilizing these classes of lipids [34,35]. Hexane, the most widely used solvent for lipid removal, although highly flammable, is an effective extraction solvent for nonpolar/neutral lipids [34,36], while reducing matrix effects [30,37]. In a study by Zhang et al. [8], acetonitrile combined with hexane liquid-liquid extraction was used to remove proteins and fats from pork and milk. Similarly, in another study by Zhang et al. [38], n-hexane was employed to remove fat from milk, eggs, and porcine muscle matrices with highly satisfactory outcomes. Thus, an acetonitrile-hexane saturated with an acetonitrile liquid-liquid extraction step was added to the pork matrix to eliminate large amounts of lipids [34]. Satisfactory recovery was achieved in this process.

O-aminobenzoic acid was not detected in the liver due to its inability to be linear in the liver matrix, thereby precluding quantification. Additionally, no blank liver matrix was identified. This analysis may also

consider the interference from endogenous substances, resulting from the complex composition of the liver matrix. The matrix is composed of a multitude of endogenous substances, including proteins, enzymes, and metabolites. These substances may interact with o-aminobenzoic acid in varied ways. For instance, o-aminobenzoic acid may bind with liver proteins through hydrogen bonding, hydrophobic interaction, etc. Such interactions modify its behavior in the detection system, and hence, a linear relationship cannot be obtained between the concentration and the response signals. In other cases, metabolic transformations may also occur. It should be noted that the liver is rich in a variety of metabolic enzymes, such as the cytochrome P450 enzyme system. Following the entry of o-aminobenzoic acid into the liver matrix, it may undergo rapid catalysis with the aid of these enzymes. This results in metabolic transformations, such as the oxidation of the amino group and the hydroxylation of the benzene ring. Consequently, the metabolic processes alter the detected concentration of o-aminobenzoic acid, preventing the linear relationship between the concentration and the detection signal.

Table 4

Validation parameters and MRLs reference regulation of target veterinary drugs in beef, pork, sheep, horse, chicken, prawn, fish, liver, milk, fat.

No.	Analyte	Matrix	Matrix-matched calibration curve	Correlation coefficient(R ²)	Linear range/(μ g/L)	LOD/(μ g/kg)	LOQ/(μ g/kg)	MRLs reference regulation
1	Acetanilide	Beef	$Y = 12,924.03X - 6334.18$	0.999	0.5~50	0.07	0.22	Pig, Horse, Cattle, Sheep, Goat Muscle: 0.01 mg/kg; Milk: 0.01 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general (Korean) https://residue.foodsafetykorea.go.kr/vd/mrl
		Pork	$Y = 29,574.1X - 7721.23$	0.992	1~50	0.09	0.30	
		Sheep	$Y = 55,623.2X + 16,368.8$	0.999	0.5~50	0.04	0.12	
		Horse	$Y = 52,837.3X + 22,144.1$	0.997	0.5~50	0.03	0.12	
		Chicken	$Y = 24,353.8X + 75,795.2$	0.995	0.5~50	0.04	0.15	
		Prawn	$Y = 7768.35X - 2955.34$	0.993	1~50	0.13	0.42	
		Fish	$Y = 19,834.3X + 1309.34$	0.998	1~50	0.11	0.37	
		Liver	$Y = 31,928.4X + 14,333.1$	0.994	0.5~50	0.06	0.20	
		Milk	$Y = 58,934.3X + 43,451.2$	0.999	0.5~50	0.03	0.10	
		Fat	$Y = 68,368.8X + 50,030.5$	0.999	0.5~50	0.02	0.06	
2	o-Aminobenzoic acid	Beef	$Y = 6983.77X + 73,455.6$	0.993	0.5~50	0.03	0.09	Pig Muscle: 0.01 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general (Korean) https://residue.foodsafetykorea.go.kr/vd/mrl
		Pork	$Y = 7938.26X + 25,469.2$	0.993	1~50	0.12	0.41	
		Sheep	$Y = 11,356.3X + 234,554$	0.999	2~50	0.18	0.62	
		Horse	$Y = 11,981.1X + 96,009.3$	0.996	0.5~50	0.07	0.23	
		Chicken	$Y = 6425.02X + 18,739.3$	0.993	1~50	0.16	0.54	
		Prawn	$Y = 3483.35X + 9345.87$	0.994	1~50	0.16	0.53	
		Fish	$Y = 7309.32X + 39,285.6$	0.992	1~50	0.10	0.33	
		Liver	/	/	/	/	/	
		Milk	$Y = 11,021.5X + 198,743$	0.995	0.5~50	0.06	0.20	
		Fat	$Y = 13,364.9X + 19,809$	0.998	1~50	0.09	0.31	
3	Pentetrazol	Beef	$Y = 8234.43X + 8782.33$	0.999	1~50	0.10	0.35	Pig, Horse, Cattle Muscle: 0.01 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general (Korean) https://residue.foodsafetykorea.go.kr/vd/mrl
		Pork	$Y = 6543.33X - 532.384$	0.997	0.5~50	0.06	0.22	
		Sheep	$Y = 16,748.8X + 9833.29$	0.999	2~50	0.25	0.84	
		Horse	$Y = 14,537.4X + 18,384.3$	0.994	2~50	0.23	0.75	
		Chicken	$Y = 7473.45X + 38,921.3$	0.994	1~50	0.09	0.31	
		Prawn	$Y = 3974.87X + 12,242.5$	0.997	2~50	0.21	0.68	
		Fish	$Y = 5476.86X + 6393.44$	0.995	1~50	0.13	0.45	
		Liver	$Y = 7783.22X + 4894.36$	0.999	3~50	0.46	1.53	
		Milk	$Y = 18,739.2X - 3338.35$	0.999	2~50	0.25	0.83	
		Fat	$Y = 23,293.7X + 2344.54$	0.999	1~50	0.13	0.43	
4	Methyl ephedrine hydrochloride	Beef	$Y = 5283.59X + 9123.67$	0.998	0.5~50	0.04	0.12	Pig, Horse, Cattle, Sheep, Goat Muscle: 0.01mg/kg; Milk: 0.005 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general (Korean) https://residue.foodsafetykorea.go.kr/vd/mrl
		Pork	$Y = 3928.5X + 2896.88$	0.996	0.5~50	0.05	0.18	
		Sheep	$Y = 18,274.8X - 2122.61$	0.999	0.5~50	0.05	0.17	
		Horse	$Y = 8763.08X + 672.34$	0.997	0.5~50	0.05	0.18	
		Chicken	$Y = 6473.23X - 1673.57$	0.996	2~50	0.18	0.60	
		Prawn	$Y = 7364.23X + 11,245.2$	0.993	0.5~50	0.04	0.14	
		Fish	$Y = 7683.7X - 2456.78$	0.999	0.5~50	0.04	0.14	
		Liver	$Y = 3789.71X - 1778.63$	0.993	2~50	0.29	0.97	
		Milk	$Y = 18,739.3X - 8123.89$	0.999	0.5~50	0.03	0.09	
		Fat	$Y = 39,893.5X - 2332.74$	0.999	0.5~50	0.02	0.05	
5	Phenacetin	Beef	$Y = 12,783.1X - 14,238.1$	0.998	0.5~50	0.03	0.10	Chicken, Pig, Horse, Deer, Cattle, Sheep, Goat, Duck, Turkey, Rabbit Muscle: 0.01mg/kg; Poultry Eggs: 0.01mg/kg; Milk: 0.01mg/kg; Source: Applicability of Veterinary Drug MRLs for food in general (Korean) https://residue.foodsafetykorea.go.kr/vd/mrl
		Pork	$Y = 17,849.4X - 8877.55$	0.997	0.5~50	0.05	0.16	
		Sheep	$Y = 89,837.1X - 2337.6$	0.999	0.5~50	0.02	0.06	
		Horse	$Y = 79,803.2X - 9762.95$	0.995	0.5~50	0.02	0.05	
		Chicken	$Y = 14,757.32X - 18,780.3$	0.995	0.5~50	0.05	0.17	
		Prawn	$Y = 9793.84X - 387.554$	0.994	5~50	0.88	2.94	
		Fish	$Y = 15,243.42X - 28,767.6$	0.992	5~50	0.85	2.84	
		Liver	$Y = 49830X + 10,983.1$	0.997	0.5~50	0.03	0.10	
		Milk	$Y = 68,973.7X - 7674.56$	0.999	0.5~50	0.02	0.08	

(continued on next page)

Table 4 (continued)

No.	Analyte	Matrix	Matrix-matched calibration curve	Correlation coefficient(R ²)	Linear range/(μ g/L)	LOD/(μ g/kg)	LOQ/(μ g/kg)	MRLs reference regulation
6	Antipyrine	Fat	$Y = 118238X - 4333.14$	0.999	0.5~50	0.01	0.04	Pig, Horse, Cattle, Sheep, Goat Muscle: 0.01 mg/kg; Milk: 0.01 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general (Korean) https://residue.foodsafetykorea.go.kr/vd/mrl
		Beef	$Y = 7832.43X + 1976.64$	0.997	0.5~50	0.08	0.26	
		Pork	$Y = 33,243.9X + 56,752.9$	0.995	0.5~50	0.02	0.06	
		Sheep	$Y = 46,453.3X - 8909.54$	0.994	0.5~50	0.02	0.06	
		Horse	$Y = 44,354.1X + 7808.32$	0.999	0.5~50	0.02	0.05	
		Chicken	$Y = 31,232.5X + 8797.73$	0.998	0.5~50	0.04	0.12	
		Prawn	$Y = 7342.65X + 2309.43$	0.997	0.5~50	0.04	0.14	
		Fish	$Y = 13,245.7X - 7909.6$	0.997	0.5~50	0.03	0.12	
		Liver	$Y = 27,242.2X + 5098.23$	0.993	0.5~50	0.05	0.16	
		Milk	$Y = 32,343.3X + 2654.89$	0.998	0.5~50	0.02	0.07	
7	Guaifenesin	Fat	$Y = 54,675.5X - 4984.45$	0.998	0.5~50	0.02	0.05	Chicken, Pig, Horse, Cattle Muscle: 0.01 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general (Korean) https://residue.foodsafetykorea.go.kr/vd/mrl
		Beef	$Y = 5865.13X - 209.453$	0.993	3~50	0.52	1.74	
		Pork	$Y = 1453.32X - 988.032$	0.996	3~50	0.43	1.45	
		Sheep	$Y = 7435.22X + 2988.65$	0.994	3~50	0.31	1.03	
		Horse	$Y = 6675.71X + 3767.77$	0.999	3~50	0.40	1.33	
		Chicken	$Y = 2149.42X + 9943.3$	0.999	5~50	0.74	2.45	
		Prawn	$Y = 7345.96X - 715.209$	0.999	5~50	0.63	2.10	
		Fish	$Y = 1764.43X - 2123.69$	0.994	5~50	0.68	2.27	
		Liver	$Y = 4324.6X - 885.232$	0.999	5~50	0.83	2.78	
		Milk	$Y = 7453.32X - 402.143$	0.999	2~50	0.29	0.96	
8	Diethylcarbamazine	Fat	$Y = 11,234.8X - 2541.06$	0.999	2~50	0.22	0.72	Horse, Cattle, Sheep, Goat Muscle: 0.01 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general (Korean) https://residue.foodsafetykorea.go.kr/vd/mr
		Beef	$Y = 25,432.8X + 178,633$	0.999	0.5~50	0.01	0.03	
		Pork	$Y = 37,643.2X + 3789.45$	0.995	0.5~50	0.01	0.03	
		Sheep	$Y = 45,643.4X - 5981.2$	0.999	0.5~50	0.03	0.10	
		Horse	$Y = 59,869.9X + 17,867.5$	0.998	0.5~50	0.02	0.07	
		Chicken	$Y = 65,436.5X + 139,733$	0.999	0.5~50	0.03	0.11	
		Prawn	$Y = 75,355.7X + 208,945$	0.995	0.5~50	0.01	0.03	
		Fish	$Y = 37,549.1X + 8554.53$	0.993	0.5~50	0.02	0.06	
		Liver	$Y = 65,413.8X - 5234.61$	0.993	0.5~50	0.02	0.05	
		Milk	$Y = 134636X - 48,891.1$	0.999	0.5~50	0.01	0.02	
9	3,6-diamino-10-methylacridinium chloride	Fat	$Y = 223,411.3X + 16,230.9$	0.998	0.5~50	0.01	0.02	Horse, Cattle, Sheep Muscle: 0.01 mg/kg; Milk: 0.01 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general (Korean) https://residue.foodsafetykorea.go.kr/vd/mrl
		Beef	$Y = 4342.87X + 1323.21$	0.992	1~50	0.17	0.50	
		Pork	$Y = 7324.43X - 402.334$	0.992	0.5~50	0.01	0.04	
		Sheep	$Y = 19,876.7X - 2543.87$	0.999	0.5~50	0.02	0.07	
		Horse	$Y = 12,321.3X - 343.724$	0.997	0.5~50	0.03	0.12	
		Chicken	$Y = 7342.89X - 7342.72$	0.995	1~50	0.10	0.33	
		Prawn	$Y = 5342.23X + 7897.66$	0.998	1~50	0.10	0.35	
		Fish	$Y = 6342.95X + 3972.34$	0.993	0.5~50	0.07	0.25	
		Liver	$Y = 4843.99X + 5698.25$	0.995	0.5~50	0.08	0.27	
		Milk	$Y = 15,462.9X - 6988.53$	0.997	0.5~50	0.05	0.17	
10	3,6-Acridinediamine	Fat	$Y = 23,343.2X - 4874.78$	0.998	0.5~50	0.01	0.04	Pig Muscle: 0.01 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general(Korean) https://residue.foodsafetykorea.go.kr/vd/mrl
		Beef	$Y = 5231.24X + 387.564$	0.999	1~50	0.13	0.43	
		Pork	$Y = 7651.32X + 3987.21$	0.998	0.5~50	0.03	0.10	
		Sheep	$Y = 24,533.6X - 1987.72$	0.999	0.5~50	0.03	0.09	
		Horse	$Y = 24,324.6X - 711.89$	0.994	0.5~50	0.02	0.08	
		Chicken	$Y = 15,343.2X + 38,973.1$	0.998	1~50	0.11	0.37	
		Prawn	$Y = 6451.16X + 4897.34$	0.999	1~50	0.11	0.35	
		Fish	$Y = 7923.67X - 5782.13$	0.998	1~50	0.08	0.26	
		Liver	$Y = 9874.46X + 3768.71$	0.998	1~50	0.15	0.49	
		Milk	$Y = 26,548.6X - 8972.19$	0.999	0.5~50	0.03	0.09	
11	Metomidate	Fat	$Y = 36,433.2X + 9553.8$	0.999	0.5~50	0.01	0.05	Pig Muscle: 0.01 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general(Korean) https://residue.foodsafetykorea.go.kr/vd/mrl
		Beef	$Y = 23,255.2X - 49,807.9$	0.998	0.5~50	0.02	0.06	
		Pork	$Y = 48,984.3X - 46,578$	0.994	0.5~50	0.01	0.03	
		Sheep	$Y = 229,381.2X - 20,391.1$	0.999	0.5~50	0.01	0.02	
		Horse	$Y = 218,307.2X - 11,135.2$	0.998	0.5~50	0.01	0.02	
		Chicken	$Y = 38,324.1X - 19,875.3$	0.991	0.5~50	0.01	0.04	
		Prawn	$Y = 24,342.2X + 28,795.9$	0.999	0.5~50	0.08	0.26	
		Fish	$Y = 21,732.7X - 34,098$	0.996	0.5~50	0.01	0.04	

(continued on next page)

Table 4 (continued)

No.	Analyte	Matrix	Matrix-matched calibration curve	Correlation coefficient(R ²)	Linear range/(μg/L)	LOD/(μg/kg)	LOQ/(μg/kg)	MRLs reference regulation
12	Ethacridine lactate monohydrate	Liver	Y = 132458X+2223.23	0.997	0.5~50	0.01	0.02	Pig, Horse, Cattle Muscle: 0.01 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general (Korean) https://residue.foodsafetykorea.go.kr/vd/mrl
		Milk	Y = 167,434.3X-57,686.3	0.999	0.5~50	0.01	0.03	
		Fat	Y = 255332X-17,817.2	0.998	0.5~50	0.01	0.03	
		Beef	Y = 26,732.2X+737.092	0.999	0.5~50	0.01	0.05	
		Pork	Y = 40,342.6X-6094.42	0.996	0.5~50	0.01	0.03	
		Sheep	Y = 130293X-19,878.6	0.999	0.5~50	0.01	0.03	
		Horse	Y = 149303X-2898.1	0.998	0.5~50	0.01	0.02	
		Chicken	Y = 33,442.6X-26,664.3	0.995	0.5~50	0.03	0.10	
		Prawn	Y = 23,244.9X+28,904	0.993	0.5~50	0.01	0.03	
		Fish	Y = 36,553.3X-56,375.5	0.994	0.5~50	0.02	0.08	
		Liver	Y = 49,736.8X+25,986.4	0.996	0.5~50	0.01	0.03	
		Milk	Y = 149280X-48,277.8	0.994	0.5~50	0.01	0.02	
13	Tripelethamine hydrochloride	Fat	Y = 198,374.1X-17,866.6	0.998	0.5~50	0.01	0.02	Pig, Horse, Cattle Muscle: 0.01 mg/kg; Milk: 0.01 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general (Korean) https://residue.foodsafetykorea.go.kr/vd/mrl Edible tissues of Cattle: 0.2 mg/kg; Cattle milk: 0.02 mg/kg Source: Code of Federal Regulations Title 21(US Food and Drug Administration) www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfCFR/CFRSearch.cfm?CFRPart=556&showFR=1
		Beef	Y = 47,638.4X+3873.5	0.997	0.5~50	0.01	0.04	
		Pork	Y = 68,394.1X-4034.3	0.994	0.5~50	0.01	0.03	
		Sheep	Y = 224,221.2X-4164.3	0.999	0.5~50	0.01	0.04	
		Horse	Y = 209238X-10,989.8	0.998	0.5~50	0.02	0.06	
		Chicken	Y = 97,384.7X+189,775	0.997	0.5~50	0.06	0.21	
		Prawn	Y = 41,211.8X+49,811.2	0.994	0.5~50	0.01	0.04	
		Fish	Y = 58,927.5X-34,532.7	0.998	0.5~50	0.01	0.04	
		Liver	Y = 49,383.7X+28,934.3	0.991	0.5~50	0.02	0.06	
		Milk	Y = 183948X-54,453	0.999	0.5~50	0.01	0.05	
		Fat	Y = 283,943.2X-267.986	0.996	0.5~50	0.01	0.03	
14	Ormetoprim	Beef	Y = 10,096.23X+21,133.7	0.999	0.5~50	0.02	0.07	Muscle and skin of salmonids in natural proportions: 0.1mg/kg Source: Maximum Residue Limits (MRLs) for Veterinary Drugs in Foods (Canada) https://www.canada.ca/en/health-canada/services/drugs-health-products/veterinary-drugs/maximum-residue-limits-mrls/list-maximum-residue-limits-mrls-veterinary-drugs-foods.html
		Pork	Y = 14,533.3X-1360.45	0.998	0.5~50	0.01	0.03	
		Sheep	Y = 8453.23X+378.123	0.997	0.5~50	0.16	0.53	
		Horse	Y = 6782.43X+354.100	0.999	3~50	0.37	1.24	
		Chicken	Y = 21,537.4X+48,938.2	0.998	0.5~50	0.03	0.09	
		Prawn	Y = 14,539.3X+6549.7	0.993	0.5~50	0.03	0.10	
		Fish	Y = 16,758.8X+4343.32	0.999	0.5~50	0.03	0.12	
		Liver	Y = 5340.64X+694.122	0.992	1~50	0.15	0.50	
		Milk	Y = 20,246.1X-5193.2	0.996	1~50	0.08	0.28	
		Fat	Y = 30,322.2X-2123.46	0.999	0.5~50	0.04	0.12	
		Beef	Y = 13,452.6X+2389.37	0.999	0.5~50	0.01	0.04	
		Pork	Y = 17,654.4X-7776.68	0.998	0.5~50	0.02	0.07	
15	Sulfaethoxypyridazine	Sheep	Y = 53,420.4X-8123.56	0.997	0.5~50	0.01	0.04	Fat, kidney, Liver, Muscle of cattle; Kidney, Liver, Muscle, Skin and fat of swine: 0.1mg/kg - singly or in combination with other sulfonamides listed. Milk of cattle: 0.01mg/kg - singly or in combination with other sulfonamides listed. Source: Maximum Residue Limits (MRLs) for Veterinary Drugs in Foods (Canada) https://www.canada.ca/en/health-canada/services/drugs-health-products/veterinary-drugs/maximum-residue-limits-mrls/list-maximum-residue-limits-mrls-veterinary-drugs-foods.html Edible tissues (excluding milk): 0.1 mg/kg; Milk: Zero; Swine. Edible tissues: Zero. Source: Code of Federal Regulations Title 21(US Food and Drug Administration) www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfCFR/CFRSearch.cfm?CFRPart=556&showFR=1
		Horse	Y = 50,424.2X-3111.64	0.998	0.5~50	0.01	0.04	
		Chicken	Y = 28,764.7X+63,452	0.992	0.5~50	0.04	0.15	
		Prawn	Y = 9975.54X+21,523.1	0.996	0.5~50	0.06	0.21	
		Fish	Y = 15,439.1X-756.657	0.999	0.5~50	0.05	0.17	
		Liver	Y = 16,532.2X+23,576.8	0.998	0.5~50	0.03	0.08	
		Milk	Y = 53,424.5X-25,779.9	0.999	0.5~50	0.01	0.04	
		Fat	Y = 79,342.5X+27,774.5	0.998	0.5~50	0.01	0.02	
		Beef	Y = 24,938.5X+21,187.3	0.997	0.5~50	0.01	0.04	
		Pork	Y = 38,971.4X+5435.2	0.998	0.5~50	0.03	0.11	
		Sheep	Y = 88,394.1X-15,647.2	0.998	0.5~50	0.01	0.03	
		Horse	Y = 76,473.2X-4988.52	0.998	0.5~50	0.01	0.03	
16	Yohimbine hydrochloride	Chicken	Y = 78,326.9X+62,222	0.999	0.5~50	0.04	0.15	Deer Muscle: 0.01 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general(Korean) https://residue.foodsafetykorea.go.kr/vd/mrl
		Prawn	Y = 21,122.4X+112,435	0.998	0.5~50	0.01	0.03	
		Fish	Y = 33,242.4X-1879.23	0.993	0.5~50	0.05	0.15	

(continued on next page)

Table 4 (continued)

No.	Analyte	Matrix	Matrix-matched calibration curve	Correlation coefficient(R ²)	Linear range/(μg/L)	LOD/(μg/kg)	LOQ/(μg/kg)	MRLs reference regulation
17	Buquinolate	Liver	Y = 21,113.6X+14,561.3	0.993	0.5~50	0.02	0.06	Kidney, Liver, Skin and fat of chickens: 0.4 mg/kg Muscle of chickens: 0.1 mg/kg Source: Maximum Residue Limits (MRLs) for Veterinary Drugs in Foods (Canada) https://www.canada.ca/en/health-canada/services/drugs-health-products/veterinary-drugs/maximum-residue-limits-mrls/list-maximum-residue-limits-mrls-veterinary-drugs-foods.html
		Milk	Y = 98,337.6X-49,876.9	0.995	0.5~50	0.01	0.03	
		Fat	Y = 133459X+2777.25	0.999	0.5~50	0.01	0.03	
		Beef	Y = 13,242.4X-189.335	0.995	0.5~50	0.07	0.24	
		Pork	Y = 63,243.6X+18,971.6	0.992	0.5~50	0.06	0.20	
		Sheep	Y = 155432X-33,457.1	0.999	0.5~50	0.01	0.03	
		Horse	Y = 85,334.7X-9878.41	0.996	0.5~50	0.01	0.02	
		Chicken	Y = 41,234.2X+18,666.8	0.997	10~50	1.26	4.20	
		Prawn	Y = 34,345.1X+15,779.3	0.999	3~50	0.50	1.68	
		Fish	Y = 64,242.1X-1223.12	0.997	0.5~50	0.06	0.22	
18	Dehydrocholic acid	Liver	Y = 58,758.7X+7534.88	0.996	0.5~50	0.01	0.03	Fish: 0.01 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general(Korean) https://residue.foodsafetykorea.go.kr/vd/mrl
		Milk	Y = 104246X-31,231.2	0.999	0.5~50	0.01	0.02	
		Fat	Y = 135334X+564.735	0.999	0.5~50	0.01	0.03	
		Beef	Y = 3346.14X-278.712	0.995	3~50	0.43	1.43	
		Pork	Y = 1768.31X-454.567	0.996	10~50	1.29	4.31	
		Sheep	Y = 4340.11X-945.762	0.999	2~50	0.32	1.07	
		Horse	Y = 3976.34X-1564.67	0.998	5~50	0.69	2.30	
		Chicken	Y = 2568.37X+3432.08	0.993	5~50	0.69	2.31	
		Prawn	Y = 5785.37X+567.098	0.997	10~50	0.88	2.94	
		Fish	Y = 2345.15X-218.225	0.994	10~50	0.83	2.78	
19	Loperamide hydrochloride	Liver	Y = 2123.21X-543.334	0.999	10~50	0.91	3.05	Chicken, Pig, Cattle Muscle: 0.01 mg/kg Source: Applicability of Veterinary Drug MRLs for food in general(Korean) https://residue.foodsafetykorea.go.kr/vd/mrl
		Milk	Y = 2531.11X-521.567	0.991	3~50	0.37	1.25	
		Fat	Y = 4345.45X-967.204	0.997	3~50	0.37	1.23	
		Beef	Y = 50,962.1X-2112.34	0.999	0.5~50	0.01	0.04	
		Pork	Y = 84,532.2X-34,347.3	0.993	0.5~50	0.02	0.06	
		Sheep	Y = 414243X-64,341.5	0.999	0.5~50	0.01	0.03	
		Horse	Y = 371461X-1435.7	0.998	0.5~50	0.01	0.04	
		Chicken	Y = 155322X+245,609	0.998	0.5~50	0.01	0.04	
		Prawn	Y = 55,740.1X+187,675	0.998	0.5~50	0.02	0.08	
		Fish	Y = 49,234.1X-69,821.9	0.994	0.5~50	0.02	0.07	
20	Ciclesonide	Liver	Y = 193241X+48,977.5	0.997	0.5~50	0.01	0.03	Muscle, Liver, Kidney of Equidae: 0.6 μg/kg; Fat of Equidae: 4 μg/kg; Not for use in animals from which milk is produced for human consumption. Source: Regulation (EU) No 43/2020; Regulation (EC) No 470/2009
		Milk	Y = 265309X-91,123	0.993	0.5~50	0.02	0.05	
		Fat	Y = 425632X+89,856.4	0.999	0.5~50	0.02	0.07	
		Beef	Y = 4532.42X-309.819	0.993	0.5~50	0.15	0.50	
		Pork	Y = 6342.17X-235.345	0.996	0.5~50	0.15	0.60	
		Sheep	Y = 4344.23X-563.074	0.999	0.5~50	0.15	0.60	
		Horse	Y = 4453.22X-193.321	0.997	0.5~50	0.15	0.60	
		Chicken	Y = 7543.64X-112.133	0.992	0.5~50	0.14	0.47	
		Prawn	Y = 4764.21X-242.456	0.997	0.5~50	0.03	0.10	
		Fish	Y = 3985.23X-135.411	0.996	0.5~50	0.15	0.50	
		Liver	/	/	/	/	/	
		Milk	Y = 8982.1X-10,788.2	0.993	0.5~50	0.16	0.52	
		Fat	Y = 6563.31X+1358.46	0.998	0.5~50	0.16	0.53	

Table 5

Peak areas, noise signals and spiked concentrations for some typical substances.

Matrix	Analyte	Peak area	Noise signal	Spiked concentration(μg/kg)	Test concentration(μg/kg)
Chicken	Sulfaethoxypyridazine	186,501.375	610.77	5	4.49
	Buquinolate	6550.171	11.937	5	5.32
Beef	Sulfaethoxypyridazine	171,544.5	2386.815	10	8.63
	Buquinolate	20,621.824	426.114	10	8.38

3.4.2. Optimization of centrifugal temperature

The samples were subjected to centrifugation at room temperature (20 ± 2 °C), at 4 °C, and at −4 °C, respectively, and the results are shown in Fig. 3. In fish (grass carp), for instance, there was a negligible difference in the recoveries of the target veterinary residues at 4 °C and at −4 °C. The recoveries of o-Aminobenzoic acid, antipyrine, and loperamide at −4 °C were found to be slightly higher than those at 4 °C, although the deviation did not exceed 2%. However, the recovery efficiency of a majority of the compounds centrifuged at room temperature exhibited a slight decline in comparison to the latter two, with

ormetoprim demonstrating notable decreases (< 9%). Centrifuges generate a considerable amount of heat when operated at high speeds, and proteins or fats present in centrifuged samples may undergo qualitative changes at elevated temperatures. According to Schmid et al. [39], the difference in temperature demonstrates a significant effect on fats, and their composition changes over time. Furthermore, Escorsim et al. [34] also states that as temperature rises, fat extraction becomes more facile, which will impact the outcome of the experiments and introduce interference with the target to be measured. Given that there was minimal distinction between the recoveries obtained at the




Fig. 7. Total ion chromatograms (TIC) and extracted ion chromatography (EIC) of sulfaethoxypyridazine (spiked conc. 5 µg/kg) in chicken(a, b, c) and (spiked conc. 10 µg/kg) in beef (d, e, f).




Fig. 7. (continued).

temperatures of 4 °C and −4 °C, the final parameter decided for the experiment was cryo-centrifugation at 4 °C.

3.5. Optimization of the purification step

3.5.1. Selection of purification powder and C18 dosage

The composition of animal food matrices is intricate, comprising a multitude of fats and phospholipids that are readily extracted by acetonitrile. However, direct injection into the column may result in adsorption and subsequent difficulty in elution. Long-term utilisation of the system may result in increased stress, necessitating the removal of impurities such as fats and phospholipids prior to mass spectrometry. As a typical dispersive solid-phase extraction (DSPE) technique, QuEChERS employs the interaction between the adsorbent filler and the interfering matrix to achieve the purpose of decontamination and clean-up. This method exhibits the advantages of simple operation, high recovery, accurate detection results, and a wide range of applications [10]. The principal purification materials employed in the QuEChERS technique include N-(n-propyl)ethylenediamine (PSA), octadecyl-bonded silica gel (C18), graphitised carbon black (GCB), anhydrous magnesium sulphate (MgSO_4), etc. PSA is employed for the removal of organic acids, pigments, and metal ions [40]. On the other hand, C18 has a greater capacity to remove non-polar substances and can effectively remove long-chain lipids, such as fats and fat-soluble pigments [41]. GCB is primarily utilized to adsorb substances such as pigments [42]. A series of experiments were conducted with the objective of comparing the efficacy of different purified powder formulations at varying ratios: 1) 80 mg C18+1000 mg MgSO_4 +80 mg PSA; 2) 80 mg C18+1000 mg MgSO_4 +40 mg GCB; 3) 80 mg C18+1000 mg MgSO_4 +80 mg PSA+40 mg GCB; and 4) 80 mg C18+1000 mg MgSO_4 , as shown in Fig. 4. The outcomes demonstrated that the inclusion of PSA enhanced the recovery rate of certain substances to a certain degree. Nevertheless, its effect was less pronounced than that of the group 4 with respect to the overall

recovery rate of substances. The incorporation of GCB notably decreased the recuperation of several substances, with a distinct impact on 3,6-diamino-10-methylacridinium and 3,6-acridinediamine. Consequently, C18 and MgSO_4 were selected as purification reagents, and their dosages were optimized. A mixture with a concentration level of 10 ng/ml was added to beef, pork, chicken, prawn, and fish under the same conditions to investigate the purification effect of the target species of veterinary residues at 40 mg, 60 mg, 80 mg, 100 mg, and 120 mg C18, respectively. Six replicates were conducted for each matrix, and the mean recoveries of the ten matrices are presented in Fig. 5. The results demonstrated that a majority of the targets exhibited superior recoveries and purification in ten representative matrices when C18 was 100 mg. Conversely, the recoveries were significantly lower when the amount of C18 powder used was 40 mg or less. The beef matrix (Fig. 5a) was used as an example to illustrate the results. The data indicated that pentetrazol, methyl ephedrine hydrochloride, and phenacetin exhibited slightly higher recoveries (not exceeding 4%) at 80 mg C18 than at 100 mg C18. Loperamide demonstrated the highest recovery at 60 mg C18, while ciclesonide exhibited the highest recovery at 120 mg C18. When C18 was set at 40 mg, the majority of substances exhibited reduced recoveries, with the exception of buquinolate and antipyrine.

3.5.2. Dosage of mgso_4

The presence of water reduces the rate of concentration. Extracted by the QuEChERS method, salt is usually added at the time of extraction in order to reduce the solubility of the veterinary drug in the aqueous phase by means of salting out. Therefore, the content of the veterinary drug in the organic phase is increased, facilitating the separation of the aqueous phase from the organic phase and reducing the amount of water in the organic phase [33]. The most commonly used salts are MgSO_4 , Na_2SO_4 and NaCl. Na_2SO_4 is particularly challenging as it can easily release a large amount of heat when exposed to water. However, veterinary residues should be kept at a relatively low temperature during the




Fig. 8. TIC and EIC of buquinolate (spiked conc. 5 µg/kg) in chicken(a, b, c) and (spiked conc. 10 µg/kg) in beef (d, e, f)




Fig. 8. (continued).

extraction process in order to minimise fat leaching [34]. Furthermore, sodium chloride may remove some water-soluble compounds. In this regard, Schenck et al. [43] have reported that MgSO_4 is a more effective drying agent than Na_2SO_4 with respect to the removal of residual water from organic solvents. Following a comparison of the various dehydrating salts, MgSO_4 was finally selected as the dehydrating salt in this study, and its content was optimised (in Fig. 6). Under the conditions of 1100 mg of MgSO_4 , 85% acetonitrile water achieved the maximum extraction efficiency as the extraction solvent, while also removing water rapidly and effectively reducing the nitrogen blowing time. Consequently, the final dosage for the purification step was determined to be 100 mg of C18 and 1100 mg of MgSO_4 .

3.6. Evaluation of the matrix effects

Following the extraction and purification of the samples, it is possible that some common extracts may affect the ionisation efficiency of the target veterinary residues. This phenomenon may result in the enhancement or inhibition of the signal response of the target substances, which is known as the matrix effect (ME). During the electrospray process, compounds may be affected by ionisation and the transfer of impurities, coeluting analytes, or degradation products from the formation of charged droplets from the sample solution ejected by the capillary tip to the entry of charged ions into the cone pore [44,45]. Matrix compounds can modify the physical properties of the droplet, including surface tension. A competition between the matrix and the analyte for the limited elemental charge on the droplet may occur. This phenomenon can interfere with the detection or recovery of the target substance [45]. This interference manifests as reduced recovery or response. The matrix enhancement or inhibition effect may also cause false-positive or false-negative results.

Consequently, the majority of researchers endeavour to employ a multitude of techniques to minimise or negate the influence of the matrix. These methods include optimisation of sample preparation [46,47],

adjustment of chromatographic parameters [48,49], utilisation of internal standards [50], setting of flow splitting [51] or the necessity to employ standard additions [44,52,53]. Therefore, this study conducted a preliminary investigation into the matrix effect of the method. The negative samples were extracted and purified in accordance with Section 2.3, resulting in the production of the blank matrix extract. The blank matrix extract and methanol were employed in the preparation of mixed standard solutions, which were used to generate matrix-matched calibration solutions and solvent standard solutions. These solutions were analyzed using the same method. The MEs were evaluated by using the following equation:

$$ME = \left(\frac{A}{B} - 1 \right) \times 100\%$$

where A is the slope of the matrix standard curve and B stands for the slope of the pure solvent standard curve.

The matrix effect varies for different kinds of substances, with $ME < 0$ representing matrix inhibition and vice versa for matrix enhancement. A larger absolute value indicates a stronger matrix effect. When $|ME| < 20\%$, a weak matrix effect is indicated, which suggests that the matrix effect is not significant; when $|ME| < 50\%$, a moderate matrix effect is indicated; and when $|ME| > 50\%$, a strong matrix effect is implied, which is considered to be significant [54]. Typical matrices employed in this study included poultry (chicken), livestock (beef, pork, sheep, horse), aquatic (Grass carp, prawn), liver, fat, and milk matrices. These matrices were selected to evaluate the matrix effects of the test targets across different product categories. As illustrated in Table 3, the targets exhibited varying degrees of matrix enhancement or matrix inhibition in all ten matrices. The target veterinary residues exhibited significant matrix inhibition in the majority of matrices, with the majority also exhibiting moderate or strong matrix effects. For instance, acetanilide, pentetrazol, methyl ephedrine hydrochloride, phenacetin, antipyrine, and diethylcarbamazine in beef; metomidate, tripeleminamine,

Table 6
Mean recoveries and RSD of the targeted veterinary residues at different spiked levels, *n* = 6.

No.	Compound	Spiked level (µg/kg)	Beef			Pork			Sheep			Horse			Chicken			Prawn			Fish			Liver			Milk			Fat								
			Rec over y(%)	RSD(%)		Rec over y(%)	RSD(%)		Rec over y(%)	RSD(%)		Rec over y(%)	RSD(%)		Rec over y(%)	RSD(%)		Rec over y(%)	RSD(%)		Rec over y(%)	RSD(%)		Rec over y(%)	RSD(%)		Rec over y(%)	RSD(%)		Rec over y(%)	RSD(%)							
				Inter day	Intra day		Inter day	Intra day		Inter day	Intra day		Inter day	Intra day		Inter day	Intra day		Inter day	Intra day		Inter day	Intra day		Inter day	Intra day		Inter day	Intra day		Inter day	Intra day	Inter day	Intra day	Inter day	Intra day	Inter day	Intra day
1	Acetanilide	5	89.3	8.5	10.6	69.6	2.3	4.6	92.1	7.5	9.9	109.2	9.8	10.3	86.1	9.4	10.5	71.2	9.6	10.5	92.1	8.3	9.5	77.3	10.3	14.6	75.2	6.9	9.8	95.7	13.0	16.5						
		10	95.6	2.5	4.1	70.4	5.3	7.6	80.2	10.5	13.6	86.6	2.3	5.3	83.4	4.3	6.7	70.8	5.6	7.6	92.9	3.9	3.5	70.9	8.4	10.5	76.3	1.0	2.3	114.5	3.5	5.8						
		50	86.7	5.4	9.5	72.7	7.4	10.2	76.8	1.4	3.4	82.8	4.0	6.4	87.2	5.4	7.5	70.6	3.6	4.6	94.5	3.4	5.7	69.7	5.4	6.4	76.6	1.3	6.5	66.4	5.8	6.8						
2	o-Aminobenzoic acid	5	84.0	3.5	7.6	62.3	2.0	5.7	66.8	3.4	4.5	63.9	1.3	2.3	77.7	6.5	8.5	72.4	7.5	9.3	67.0	3.5	6.7	/	/	/	63.5	1.2	2.3	107.9	8.4	9.5						
		10	78.9	9.6	11.2	113.0	8.6	11.5	67.6	3.8	5.8	69.5	5.9	7.9	105.1	7.6	9.3	65.4	2.5	3.5	65.4	4.6	5.4	/	/	/	67.7	4.5	5.7	112.4	4.2	5.3						
		50	113.2	2.8	5.6	116.2	1.6	1.8	65.1	2.9	3.4	67.3	3.6	4.5	103.2	5.4	8.6	103.2	7.6	9.7	70.2	5.4	6.3	/	/	/	63.7	2.5	3.3	109.7	2.5	3.5						
3	Pentetrazol	5	72.3	10.4	12.8	74.3	11.6	12.6	64.0	1.8	3.3	113.6	3.6	4.4	73.2	5.3	7.8	71.2	8.5	9.5	87.0	7.6	8.5	79.8	16.9	20.6	81.0	3.7	4.5	99.8	16.4	17.4						
		10	88.7	6.5	8.7	80.6	4.4	5.8	68.9	4.2	5.3	110.1	7.6	8.1	82.5	8.5	9.7	69.2	5.6	7.6	89.4	7.2	9.3	73.2	8.4	10.5	84.0	2.5	3.4	95.7	3.2	5.4						
		50	82.0	6.5	9.8	73.2	3.6	6.5	70.6	3.5	6.6	92.1	4.0	5.6	80.5	9.5	10.6	73.9	6.8	8.7	112.6	2.6	5.6	65.9	2.7	3.5	79.3	3.6	6.5	65.7	2.2	3.5						
4	Methyl ephedrine hydrochloride	5	67.4	4.5	6.7	65.4	3.4	4.5	97.0	10.4	13.4	98.0	5.1	6.1	82.3	3.4	5.6	72.3	6.5	9.8	77.0	3.4	5.4	83.7	3.8	4.3	65.1	4.2	5.3	90.1	3.0	6.5						
		10	70.6	7.5	8.9	75.6	2.0	4.5	74.6	5.6	6.4	85.9	10.5	12.4	76.5	4.5	7.5	75.2	4.6	6.7	82.4	5.4	6.7	75.0	7.2	8.5	67.3	4.6	4.5	95.3	1.4	2.5						
		50	69.5	4.7	6.8	70.2	6.5	7.5	62.9	1.6	3.4	85.3	1.5	2.6	69.2	5.4	8.6	71.1	5.4	8.7	76.5	7.3	8.6	75.2	5.7	5.6	66.6	3.2	3.8	85.8	1.5	4.2						
5	Phenacetin	5	93.4	8.4	9.8	65.3	4.2	4.3	70.1	6.2	8.4	112.2	6.0	8.4	74.2	5.4	7.5	78.5	8.7	9.8	79.0	9.2	11.5	62.1	1.9	2.4	76.3	6.9	7.1	93.8	2.5	3.5						
		10	88.6	8.7	10.4	66.5	6.7	5.2	97.8	2.6	3.2	99.4	2.2	5.6	84.4	6.6	6.9	76.5	5.4	6.7	72.1	5.4	6.5	67.7	4.7	5.9	79.4	2.3	3.3	95.8	4.7	5.6						
		50	91.8	5.6	6.8	71.8	6.0	5.7	67.2	2.3	4.5	89.2	7.7	9.5	81.6	5.6	9.6	77.0	4.7	8.9	78.5	3.4	7.6	67.3	5.6	6.5	82.1	1.2	2.1	87.7	12.3	13.6						
6	Antipyrine	5	91.4	13.5	18.7	66.3	3.9	5.4	74.4	7.4	11.2	98.1	16.4	18.8	87.5	5.6	8.7	69.0	7.8	8.7	88.2	10.5	13.5	73.0	9.3	11.2	82.3	6.3	6.9	97.0	6.3	7.6						
		10	94.7	11.6	12.4	67.6	5.0	6.6	79.9	2.6	3.4	103.4	2.5	2.3	92.3	4.6	5.6	65.5	4.6	6.7	73.4	3.5	3.5	70.2	7.1	8.5	85.9	1.5	1.9	98.1	5.5	6.5						
		50	111.4	5.8	5.7	70.3	4.9	6.5	70.3	3.3	6.5	92.1	2.4	5.4	84.4	5.6	7.6	68.6	5.7	6.8	86.5	4.5	7.6	79.7	2.4	3.5	87.0	2.3	3.3	67.6	3.7	4.5						

(continued on next page)

Table 6 (continued)

7	Guaiifenesin	5	95.0	6.7	7.6	71.2	6.1	6.4	101.9	14.1	17.4	92.5	3.9	3.2	82.4	8.5	12.4	75.0	10.7	11.8	76.5	9.5	10.3	107.9	15.9	17.3	78.0	6.7	6.8	87.3	2.9	3.3
		10	96.6	8.9	9.8	76.7	4.9	5.7	73.1	6.5	8.4	95.6	12.7	15.2	89.6	6.4	8.6	75.4	5.7	6.7	84.3	4.5	5.4	71.3	7.4	11.9	85.3	2.5	3.5	89.6	2.5	3.2
		50	97.5	8.5	9.7	69.0	4.3	4.9	64.4	1.3	2.3	86.5	5.7	11.2	94.6	6.8	7.5	68.5	4.6	5.8	86.7	3.5	6.7	74.6	3.3	5.6	77.2	3.4	4.5	82.4	1.7	2.5
8	Diethylcarbamazine	5	81.6	6.7	8.9	68.3	3.0	5.4	86.6	3.9	5.4	90.5	7.0	9.6	67.4	5.8	6.8	72.5	7.5	8.6	73.0	4.5	7.6	65.4	4.0	4.6	63.4	1.2	2.3	86.5	2.7	3.5
		10	95.5	7.6	9.6	76.5	2.0	3.6	71.1	2.5	3.4	88.6	4.9	5.7	66.9	2.4	4.5	74.2	8.5	9.6	71.2	2.5	3.5	63.7	2.1	3.1	64.2	1.3	2.5	89.6	5.5	6.5
		50	81.0	8.7	7.8	67.4	4.2	5.8	67.4	2.4	5.4	78.7	2.1	7.3	73.4	4.2	7.6	66.4	2.6	3.8	74.3	3.2	6.5	64.6	2.4	2.5	66.4	2.5	3.5	64.7	1.7	1.9
9	3,6-diamino-10-methylacridinium chloride	5	70.3	7.5	8.9	68.6	3.7	5.6	88.4	10.5	13.5	65.2	2.5	3.8	71.3	15.6	18.7	69.5	4.6	5.9	67.2	5.3	7.6	64.1	3.2	4.5	65.1	2.4	3.6	88.8	1.8	2.5
		10	82.8	4.6	7.8	68.9	2.5	3.4	72.3	7.1	9.5	69.6	4.2	6.6	73.5	4.6	5.8	66.2	4.6	5.7	66.4	3.4	4.4	73.7	2.3	2.8	68.5	1.6	2.3	84.4	4.7	6.5
		50	70.0	7.7	9.5	70.2	4.0	6.4	65.3	2.2	3.4	66.5	1.2	2.5	64.2	2.5	3.6	64.3	2.8	3.6	67.4	2.6	5.8	88.0	3.3	5.4	67.5	2.3	3.5	80.5	13.3	15.6
10	3,6-Acridinediamine	5	72.4	8.5	10.8	67.5	3.7	7.6	64.6	3.5	3.6	76.4	6.0	7.0	75.3	6.4	7.7	66.0	4.9	5.7	71.4	7.3	10.8	64.6	3.3	3.8	68.4	5.6	6.5	88.5	1.2	1.5
		10	72.2	8.6	9.8	69.3	5.7	6.4	73.5	3.5	5.6	85.4	1.8	2.3	84.3	5.4	8.6	67.9	4.7	5.6	73.4	7.5	5.4	67.1	4.1	5.4	73.7	1.2	2.5	85.2	1.3	1.9
		50	72.7	4.5	7.6	64.5	2.5	4.5	70.6	5.0	6.6	69.2	1.8	1.9	70.5	5.5	9.6	69.2	6.5	6.7	82.7	7.5	9.8	89.8	7.2	9.5	72.1	1.0	3.5	76.6	3.7	4.5
11	Metomidate	5	83.4	11.5	14.9	69.1	2.8	5.7	74.6	8.6	9.8	108.5	7.0	8.8	77.6	7.2	10.5	64.0	2.4	4.9	75.3	11.5	17.6	62.3	1.9	2.5	86.0	5.4	6.5	96.0	1.8	1.3
		10	70.8	3.5	6.7	71.2	3.0	4.5	80.4	1.3	2.5	106.2	3.9	5.6	72.4	4.6	6.7	67.3	8.9	11.7	72.0	3.2	4.5	68.3	4.9	5.6	90.9	2.5	2.3	97.8	1.2	1.9
		50	83.3	6.6	8.7	67.4	5.6	6.5	69.8	5.4	6.5	96.1	1.7	3.5	81.5	7.4	8.7	69.9	6.4	7.8	76.3	4.5	7.6	75.8	5.4	7.5	87.2	1.3	3.2	85.3	4.8	5.6
12	Ethacridine lactate monohydrate	5	71.8	8.4	9.7	64.5	2.9	3.7	63.9	2.7	3.1	63.9	1.9	2.8	67.3	3.4	4.7	62.3	1.5	1.9	64.0	3.6	3.7	69.5	5.4	6.5	70.1	6.6	8.5	64.9	3.5	3.8
		10	68.9	6.5	8.8	69.6	8.7	9.8	67.6	1.2	3.5	62.5	2.3	5.9	63.4	1.7	3.3	68.9	2.5	3.7	68.4	3.3	3.5	64.4	2.5	3.4	79.4	3.5	5.6	66.0	2.6	3.5
		50	71.2	8.6	9.8	70.1	3.2	6.4	68.1	3.8	6.5	69.2	4.4	5.6	68.9	2.4	3.2	67.5	4.6	5.6	68.2	2.6	4.8	70.3	1.4	3.2	72.9	2.5	3.5	65.6	3.8	3.8
13	Tripeleennamine hydrochloride	5	82.3	5.7	7.6	72.3	6.7	8.3	91.1	6.9	8.9	113.4	3.8	7.4	84.4	9.4	10.5	68.7	5.9	7.6	87.8	9.4	12.5	68.9	5.2	6.6	66.2	1.7	3.5	86.6	3.3	5.6
		10	84.4	8.6	8.8	70.6	2.5	4.3	77.8	3.5	5.4	112.2	4.4	6.5	73.3	6.5	8.6	69.4	5.6	6.7	83.2	5.4	8.6	78.8	8.8	8.9	72.5	9.5	11.5	87.3	5.9	6.1
		50	80.2	3.6	5.7	73.4	7.5	8.6	67.2	2.8	2.6	99.1	2.7	3.5	89.5	8.6	9.7	71.3	5.8	7.6	84.5	6.4	8.7	75.3	5.7	6.4	70.5	3.5	6.5	80.3	7.4	8.9
14	Ormetoprim	5	79.8	7.6	8.9	81.7	12.4	15.6	83.0	15.4	19.6	112.1	5.4	6.6	70.2	6.5	7.7	75.4	6.8	8.6	88.6	11.3	17.5	66.5	3.1	4.6	82.9	9.3	10.5	86.1	1.8	2.5
		10	85.4	11.7	12.4	74.8	3.5	5.4	78.6	3.7	5.7	107.4	8.1	10.5	74.4	7.6	8.9	83.5	4.6	7.6	74.3	7.3	8.4	78.8	7.8	8.7	85.6	4.5	5.5	91.9	9.3	12.6
		50	84.4	9.6	11.2	76.9	5.6	6.8	75.2	5.2	6.6	93.3	4.9	6.4	82.1	9.4	10.2	73.4	8.6	9.4	81.5	6.4	8.7	87.2	13.5	15.6	91.2	3.5	5.6	77.5	2.7	3.5

(continued on next page)

Table 6 (continued)

1 5	Sulfaethoxypyridazine	5	78.7	7.8	8.8	65.6	3.1	4.0	74.8	11.7	12.6	68.2	5.3	5.9	91.7	6.5	8.7	73.2	4.7	5.6	77.6	5.3	8.7	65.9	3.3	4.6	81.0	8.8	9.8	85.7	6.0	7.6
		10	81.2	7.6	8.6	76.9	4.6	5.7	65.1	1.9	2.3	72.2	5.7	6.8	82.3	6.8	8.7	82.5	10.6	13.8	70.9	2.4	4.5	72.4	8.3	9.4	79.8	2.4	3.5	77.1	2.7	3.8
		50	78.9	7.5	7.8	81.0	3.3	5.4	63.1	1.3	3.2	67.9	1.5	2.5	86.4	8.7	9.6	76.4	6.5	7.6	76.2	5.4	6.9	69.3	6.3	8.9	90.6	7.7	8.5	72.4	5.6	6.5
1 6	Yohimbine hydrochloride	5	86.4	12.6	17.2	70.6	9.7	10.8	88.0	6.6	6.8	103.8	12.9	15.2	85.8	15.6	17.8	68.5	2.7	3.4	82.2	2.5	4.5	78.1	8.9	9.8	81.1	5.9	6.5	86.6	2.4	3.5
		10	96.3	5.6	7.8	70.0	4.0	3.5	74.8	3.3	4.5	106.4	4.0	6.2	93.4	4.6	6.7	69.4	5.5	7.6	74.9	5.5	6.5	75.5	10.2	12.3	74.4	3.6	3.5	88.9	3.5	6.5
		50	81.8	6.7	8.8	69.3	5.4	6.5	66.4	5.1	5.8	96.2	5.8	6.7	82.4	6.8	8.7	65.5	2.5	2.7	84.3	6.3	7.5	91.6	8.7	8.9	85.8	1.3	2.4	71.9	7.4	7.8
1 7	Buquinolate	5	73.5	7.6	11.5	65.7	4.5	5.4	88.5	5.3	6.5	106.9	9.3	10.5	87.5	14.5	16.9	75.5	8.6	8.7	86.6	11.7	16.7	84.2	6.3	6.7	66.1	3.0	4.6	86.3	13.3	16.6
		10	70.5	4.6	5.7	76.4	4.9	6.5	71.4	1.9	1.2	110.2	3.6	4.6	87.5	7.3	8.9	71.5	8.7	9.4	93.0	2.2	4.3	90.0	9.2	11.1	62.8	1.4	1.8	92.6	8.7	9.8
		50	76.7	6.8	7.6	74.3	6.3	7.8	68.1	1.8	3.5	92.8	1.8	2.3	104.4	6.9	8.7	70.5	6.5	7.6	95.4	3.2	5.4	104.1	5.4	6.8	64.5	2.4	3.5	87.4	3.6	6.5
1 8	Dehydrocholic acid	5	78.5	8.6	13.7	71.2	6.7	8.7	102.8	14.6	15.4	106.8	13.2	14.5	66.9	1.9	5.6	69.5	3.6	4.3	67.0	4.3	6.5	67.2	4.4	5.1	63.5	1.3	1.9	83.9	2.0	3.5
		10	90.0	10.5	11.2	77.0	7.6	9.7	85.7	11.6	13.9	101.4	4.6	6.8	71.5	7.8	8.3	68.6	5.7	6.5	70.4	6.5	8.7	68.6	4.4	6.5	67.7	4.3	5.5	97.9	11.5	15.9
		50	75.3	8.7	9.6	75.9	9.4	11.8	71.8	1.6	2.2	96.5	6.3	7.2	72.4	8.5	9.8	66.6	2.5	3.8	66.5	2.6	5.6	63.3	1.7	2.4	70.1	4.5	3.9	107.2	12.9	13.5
1 9	Loperamide hydrochloride	5	85.6	9.6	10.5	68.6	2.1	3.8	91.4	1.5	2.5	93.9	15.6	18.6	84.2	10.8	12.7	80.4	9.5	10.7	76.5	12.4	13.9	72.3	8.2	9.5	85.9	11.3	12.5	67.1	6.2	7.8
		10	75.0	6.7	7.6	69.6	4.7	5.0	81.2	1.0	2.5	112.3	1.6	5.4	83.4	7.6	8.7	85.3	5.6	8.6	83.5	8.5	9.1	87.9	11.4	12.5	91.0	3.5	3.6	68.3	3.1	4.6
		50	86.4	8.4	8.6	67.2	4.5	8.0	73.3	4.2	4.7	99.2	5.9	9.0	76.5	5.7	7.9	74.4	7.6	8.7	77.0	3.4	5.4	73.1	2.4	3.7	94.0	2.5	2.8	70.6	2.4	3.5
2 0	Ciclesonide	0.6	87.5	14.5	16.7	74.5	8.7	11.6	78.0	7.7	9.4	69.2	3.5	5.4	71.2	4.7	6.8	80.5	16.5	19.7	75.6	5.3	6.4	/	/	/	103.4	10.7	11.5	101.7	9.0	11.9
		1.2	74.0	7.6	8.7	79.6	5.7	6.6	69.8	3.2	6.6	83.5	11.0	15.4	76.4	6.6	8.7	93.6	7.8	6.8	74.0	7.4	8.9	/	/	/	106.6	9.8	10.8	106.0	1.7	2.6
		10	78.5	8.6	9.8	83.4	5.7	8.7	77.1	9.9	11.2	78.2	10.8	11.5	81.7	8.9	9.7	88.6	6.8	7.9	87.6	5.4	7.6	/	/	/	109.1	2.5	3.5	94.2	4.4	5.8

Table 7

Comparison of this method with other detection methods using LC-MS/MS.

NO.	Compound	Method	Pretreatment	LOQ($\mu\text{g/kg}$)	Recovery(%)	RSD (%)	Foodstuffs	References
1	Acetanilide, Antipyrine, Anthranilic acid, Diphenhydramine, Cyproheptadine, DL-methylephedrine, Phenacetin	LC-MS/MS	Acetonitrile-QuEChERS	0.9–7.1	71.7–102.4	≤ 8.6	Milk	[61]
2	Tetramisole, Diethylcarbamazine, Guaifenesin	LC/MS/MS	Acetonitrile	0.2–2	67.47–97.38	< 20	Milk, eggs, and porcine muscle	[38]
3	Flumethasone, DL-Methylephedrine, 2-Hydroxy-4,6-dimethylpyrimidine	LC-MS/MS	Acetonitrile-QuEChERS	2–10	73.62–112.70	≤ 20.33	Porcine muscle, and pasteurized cow milk	[8]
4	Acetanilide, Pentylenetetrazole, Phenacetin, Tetramethrin	LC/ESI-MS/MS	0.1 % formic acid in acetonitrile-QuEChERS	0.5–2.5	63.75–89.30	≤ 15.78	Porcine muscle, milk, and table eggs	Zhang et al. 2016c
5	64 compounds(Acetanilide, DL-Methylephedrine HCl, Phenacetin, Antipyrine, Diethylcarbamazine, Tripeleennamine, Yohimbine, Loperamide etc.)	LC/MS/MS	Acetonitrile/water (4:1, v/v)-QuEChERS	0.03–3	60.0–120.4	≤ 30.8	Flat fish, eel, shrimp	[56]
6	Scopolamine, Metoclopramide, Acriflavine, Berberine, Tripeleennamine, Diphenhydramine, Acrinol, Triamcinolone, Loperamide, and Roxithromycin	LC/MS/MS	0.1% formic acid acetonitrile	0.5–10	70.1–93.3	19.8	Pork, milk, and eggs	[68]
7	Sulfadimethoxine, and Ormetoprim	LC-MS/MS	Acetonitrile	100	97–99	< 16	Fish feed and fish fillets	[69]
8	Nequinatate, and Buquinolate	LC-MS/MS	Acetonitrile	0.001	89.5–108.6	< 20	Chicken muscle, chicken liver, chicken heart, swine muscle, swine heart, cattle muscle, sheep muscle, egg	[67]
9	115 Drugs	LC-MS/MS	1% (v/v) acetic acid acetonitrile-QuEChERS	0.73–745.15	67.3–117.9	< 19.6	Beef	[31]
10	Metoserpatate, Buquinolate, and Diclofenac	LC-MS/MS	0.1% formic acid in acetonitrile	1–2	74.06–108.65	13.67	Pork, milk, and eggs	[66]
11	Anesthetics, and Sedatives(Metomidate etc.)	LC-ESI/MSMS	Acetonitrile and 0.1 % ammonium acetate in acetonitrile	0.5–5	64.7–112.5	1.0–8.6	Flatfish, eel, shrimp	[70]
12	19 Compounds	LC-MS/MS	Acetonitrile-QuEChERS	0.02–4.31	60.6–117.7	≤ 20.6	Beef, pork, sheep, horse, chicken, prawn, fish, liver, milk, fat	This work

sulfaethoxypyridazine, yohimbine hydrochloride, and loperamide hydrochloride in prawn demonstrated strong matrix effects ($|\text{ME}| > 50\%$). On the other hand, o-aminobenzoic acid, ormetoprim, buquinolate, and ciclesonide showed a moderate matrix effect ($50\% > |\text{ME}| > 20\%$) in pork.

Given the inherent complexity of animal-derived food matrices, the matrix-matched standard curve external standard method was employed in this study to mitigate the impact of matrix effects on the determination of the actual samples and to enhance the accuracy of the method for quantification.

3.7. Method validation

A series of matrix standard working solutions of 0.5, 1.0, 2.0, 3.0, 5.0, 10.0, 20.0, and 50.0 ng/mL, prepared in accordance with Section 2.2.2, was utilized. The matrix-matched calibration curves of the 19 target veterinary residues were plotted. The horizontal coordinate x represented the concentration of the analyte and the vertical coordinate y represented the measured peak area. The results are shown in Table 4. The correlation coefficients (R^2) of the working curves in the range of 0.5–50.0 ng/mL were 0.991–0.999, indicating a good linear relationship.

It is important to note that the limits of detection (LODs) and limits of quantitation (LOQs) values varied across different matrices for the same analyte. The LODs and LOQs were determined by analysing the spiked samples with signal-to-noise(S/N) ratios of 3 and 10, respectively. The peak areas, noise signals, and spiked concentrations of some typical substances were shown in Table 5. Representative transition chromatograms of spiked standard of chicken and beef samples were shown in Figs. 7 and 8. The LOD in cattle, chicken, prawn, fish, pork, sheep,

hourse, liver, fat, and milk were 0.01–0.52, 0.01–1.26, 0.01–0.88, 0.01–0.85, 0.01–1.29, 0.01–0.32, 0.01–0.69, 0.01–0.91, 0.01–0.37, and 0.01–0.29 $\mu\text{g/kg}$ ($S/N \geq 3$) and the LOQ were 0.03–1.74, 0.04–4.20, 0.03–2.94, 0.04–2.84, 0.03–4.31, 0.02–1.07, 0.02–2.30, 0.02–3.05, 0.02–1.23, and 0.02–0.96 $\mu\text{g/kg}$ ($S/N \geq 10$), respectively (Table 4). Upon comparison, the LOQs of all the analytes in the method were found to be lower than the veterinary drug limits set by Korea, Canada, US, and the European Union. This indicates that the method meets the requirements for routine testing of veterinary drug residues in multi-matrix animal-derived foods. The specificity of the method was analyzed by using five blank samples of each matrix. None of the chromatographic peaks were detected at the retention times of the target analytes. These results indicate that there were no matrix compounds in the purified liquids of blank samples that might have produced a false positive signal.

In addition, the recoveries, and relative standard deviations (RSDs) were investigated to evaluate the application accuracy and precision for the determination of 19 veterinary residues. The method was validated in accordance with the CODEX guidelines(CAC/GL 71–2009) [55,56] and the Ministry of Food and Drug Safety's Guidelines for Analysis of Residual Animal Drugs of Korea [57]. Recoveries of 60–120% and relative standard deviations of 20–30% or lower served as the validation criteria for the determination of the recoveries and RSD. The spiked levels were determined according to the requirements of the limits for the target substances in the Applicability of Veterinary Drug MRLs for food in general [58], Maximum Residue Limits (MRLs) for Veterinary Drugs in Foods [59], Regulation (EU) No 43/2020, Regulation (EC) No 470/2009, and Code of Federal Regulations Title 21 [60]. The samples

were spiked with 5, 10, and 50 µg/kg of the analyte (0.6, 1.2, and 10 µg/kg for ciclesonide). Six replicates of three quality control samples were analyzed on the same day and on three consecutive days, respectively.

The results showed that the target animal residues could achieve satisfactory recovery rates (60.6–117.7 %) in ten matrices. The RSD for both intra-day and inter-day precision, shown in Table 6 ($n = 6$), were $\leq 14.5\%$ and $\leq 16.7\%$ in beef, $\leq 12.4\%$ and $\leq 15.6\%$ in pork, $\leq 15.4\%$ and $\leq 19.6\%$ in sheep, $\leq 16.4\%$ and $\leq 18.8\%$ in horse, $\leq 15.6\%$ and $\leq 18.7\%$ in chicken, $\leq 16.5\%$ and $\leq 19.7\%$ in prawn, $\leq 12.4\%$ and $\leq 17.6\%$ in fish, $\leq 16.9\%$ and $\leq 20.6\%$ in liver, $\leq 11.3\%$ and $\leq 12.5\%$ in milk, $\leq 16.4\%$ and $\leq 17.4\%$ in fat, respectively. From the collected data, it can be inferred that the developed UPLC-MS/MS method was precise and accurate.

The LC-ESI(+)-MS/MS method described in this study is capable of complying with the daily detection requirements for target veterinary residues, which is significantly lower than the maximum limit level compared with the corresponding standards. Nevertheless, the lower limit requirements for ciclesonide in the EU regulations indicate a necessity for further optimisation of the conditions and modes of detection for ciclesonide. Atmospheric pressure chemical ionisation (APCI) or atmospheric pressure photoionization (APPI) sources are employed primarily for the analysis of compounds with medium/non polarity. In certain instances, analytes are unable to generate ions with sufficient strength via electrospray ionisation (ESI) due to structural and polarity constraints. In such cases, the utilisation of APCI can enhance the ion yield, thereby serving as a valuable complement to ESI. Ciclesonide is a neutral compound, and it is difficult to detect it with adequate sensitivity (Yamamoto et al. 2021). Detection using LC-APCI-MS/MS (Su et al. 2011) or LC-APCI-MS/MS modes ([7], Mascher et al. 2008), which can assist ionization, may result in lower detection or quantitation lines.

3.8. Comparison with other reported methods

The principal parameters of the proposed method were contrasted with other LC-MS/MS methods reported in the literature for the detection of veterinary residues in foods of animal origin (Table 7). In a study conducted by Kim et al. [61], LC-MS/MS was employed to detect seven veterinary residues in milk. Five of these were found to be identical to those identified in the present study. Similarly, Zhang et al. [38] employed LC/ESI-MS/MS to detect four veterinary residues in porcine muscle, milk, and table eggs. Three of these were found to be identical to those identified in the present study. Nevertheless, the assay developed in this study is capable of achieving a lower limit of quantification and encompasses a broader range of matrices, including beef, chicken, and fish, which renders it more widely adaptable. Kim et al. [56] employed LC-MS/MS to identify 64 veterinary residues, including acetanilide, dl-methylephedrine HCl, phenacetin, antipyrine, diethylcarbamazine, tripelemnamine, yohimbine, loperamide, and others, in flatfish, eel, and shrimp. Nevertheless, this study was able to achieve lower quantification limits for a number of compounds, including acetanilide, dl-methylephedrine HCl, loperamide, and diethylcarbamazine. In order to comply with the regulatory limits set by Korea, Canada, the EU, and the USA, a number of different matrix assays were employed. The results demonstrated satisfactory accuracy and precision. Yohimbine is an α 2-adrenergic receptor antagonist that is most commonly used in veterinary medicine to reverse the effects of the α 2-receptor agonists, xylazine and detomidine [62]. However, the focus of the current detection methods is on botanicals, barks and dietary supplements [63–65], with limited research having been conducted on the detection of this substance in food matrices of animal origin. Nakajima et al. [66] developed an LC-ESI-MS/MS method for the determination of nequinat and buquinolone in eight matrices (chicken muscle, liver, heart etc.). In contrast, this method can detect a greater number of substances and involves a greater number of matrices. In a study published by Jung et al. [67], an LC-MS/MS method was developed for the determination of 115

veterinary drug residues in beef. Of these, 10 overlapped with those tested in this paper, but had higher LOQs (acetanilide: 1.2 µg/kg, dl-methylephedrine: 1.0 µg/kg, phenacetin: 1.3 µg/kg, antipyrine: 1.7 µg/kg, diethylcarbamazine: 1.2 µg/kg, acriflavine: 2.7 µg/kg, tripelemnamine: 1.18 µg/kg, ormetoprim: 1.6 µg/kg, yohimbine: 1.4 µg/kg, loperamide: 0.8 µg/kg) and did not adequately cover their matrices. The authors of the paper also noted that their LOQs were higher than those reported in previous studies of multi-residue analyses of beef. As per the discussion, it can be concluded that the QuEChERS combined with the UPLC-MS/MS method proposed in this work has more advantages over other methods.

3.9. Determination of veterinary drugs in real samples

The applicability of the optimized UPLC-MS/MS method was evaluated on authentic samples procured from local markets. A total of 50 batches (5 batches per matrix) of the samples were subjected to analysis and quantification using matrix-matched analytical calibration curves. The results demonstrated that none of the aforementioned substances were detected.

4. Conclusions

In this study, a QuEChERS coupled with UPLC-MS/MS method was developed for the determination of 19 veterinary drug residues in foods of animal origin. This method enabled the simultaneous detection of multiple veterinary residues in multiple matrices (beef, pork, sheep, horse, chicken, prawn, fish, liver, milk, fat). At present, there is a paucity of literature pertaining to the detection of 19 target veterinary drug residues in animal-derived foods. It should be noted that some countries (e.g. China) have yet to establish relevant limits and testing standards. In comparison to existing reports, the method proposed here is relatively simple, highly sensitive, and applicable to a variety of matrices. The method complies with the limit requirements of relevant matrices in the food safety standards of South Korea, the United States, Canada, and other countries. It has been successfully used to detect a variety of animal-derived foods in the local market, providing valuable technical assistance for cross-border food safety testing.

CRedit authorship contribution statement

Qianran Sun: Writing – review & editing, Writing – original draft. **Jun Liu:** Project administration, Funding acquisition. **Yuan Gou:** Formal analysis, Data curation. **Tieyuan Chen:** Investigation, Formal analysis. **Xiaofang Shen:** Methodology, Investigation. **Tao Wang:** Validation, Supervision. **Yongli Li:** Software, Resources. **Huizhen He:** Validation, Software. **Huidan Deng:** Visualization, Validation. **Yi Hua:** Visualization, Validation.

Declaration of competing interest

The authors declare that they have no known competing financial interests or personal relationships that could have appeared to influence the work reported in this paper.

Acknowledgments

This work was supported by the National Key Research and Development Program of China [grant number 2022YFF1100804].

Data availability

Data will be made available on request.

References

- [1] J. Chen, Z. Wei, X.-y. Cao, QuEChERS pretreatment combined with ultra-performance liquid chromatography–tandem mass spectrometry for the determination of four veterinary drug residues in marine products, *Food Anal. Methods* 12 (2019) 1055–1066.
- [2] L. Jank, M.T. Martins, J.B. Arsand, T.M.C. Motta, T.C. Feijó, T. dos Santos Castilhos, R.B. Hoff, F. Barreto, T.M. Pizzolato, Liquid chromatography–tandem mass spectrometry multiclass method for 46 antibiotics residues in milk and meat: development and validation, *Food Anal. Methods* 10 (2017) 2152–2164.
- [3] M.D. Marazuela, Determination of veterinary drug residues in foods by liquid chromatography–mass spectrometry: basic and cutting-edge applications, *Liquid Chromatography* (2017) 539–570.
- [4] R. Pratiwi, S.P. Ramadhanti, A. Amatullo, S. Megantara, L. Subra, Recent advances in the determination of veterinary drug residues in food, *Foods* 12 (18) (2023) 3422.
- [5] J.-W. Zhou, X.-M. Zou, S.-H. Song, G.-H. Chen, Quantum dots applied to methodology on detection of pesticide and veterinary drug residues, *J. Agric. Food Chem.* 66 (6) (2018) 1307–1319.
- [6] E. Van Asselt, J. Jager, L. Jansen, E. Hoek-van den Hil, I. Barbu, P. Rutgers, M. Pikkemaat, Prioritizing veterinary drug residues in animal products for risk-based monitoring, *Food Control* 151 (2023) 109782.
- [7] W.-S. Chae, S. Kim, H.-J. Lee, Development of an analytical method for the determination of dl-methylephedrine hydrochloride in porcine muscle using liquid chromatography–tandem mass spectrometry, *Korean J. Veterinary Res.* 60 (4) (2020) 209–213.
- [8] D. Zhang, J.-A. Park, S.-K. Kim, S.-H. Cho, D. Jeong, S.-M. Cho, H. Yi, J.-H. Shim, J.-S. Kim, A. Abd El-Aty, Simultaneous detection of flumethasone, dl-methylephedrine, and 2-hydroxy-4, 6-dimethylpyrimidine in porcine muscle and pasteurized cow milk using liquid chromatography coupled with triple-quadrupole mass spectrometry, *J. Chromatogr. B* 1012 (2016) 8–16.
- [9] N. Drabińska, M.A. Marcinkowska, M.N. Wiczorek, H.H. Jeleń, Application of sorbent-based extraction techniques in food analysis, *Molecules* 28 (24) (2023) 7985.
- [10] G. Islas, I.S. Ibarra, P. Hernandez, J.M. Miranda, A. Cepeda, Dispersive solid phase extraction for the analysis of veterinary drugs applied to food samples: a review, *Int. J. Anal. Chem.* 2017 (2017).
- [11] T. Yao, J. Song, C. Zhou, X. Shi, Recent progress of the applications of functionalized magnetic ionic liquids in sample pretreatment, *Sep. Purif. Technol.* 341 (2024) 126979.
- [12] S.M. Plakas, K.R. El Said, E.L. Jester, F.A. Bencsath, W.L. Hayton, Liquid chromatographic determination of acriflavine and proflavine residues in channel catfish muscle, *J. AOAC Int.* 80 (3) (1997) 486–490.
- [13] S.B. Turnipseed, J.M. Storey, I.-L. Wu, W.C. Andersen, M.R. Madson, Extended liquid chromatography high resolution mass spectrometry screening method for veterinary drug, pesticide and human pharmaceutical residues in aquaculture fish, *Food Addit. Contam. Part A* 36 (10) (2019) 1501–1514.
- [14] R. Bandu, H.S. Ahn, J.W. Lee, Y.W. Kim, S.H. Choi, H.J. Kim, K.P. Kim, Liquid chromatography electrospray ionization tandem mass spectrometry (LC/ESI-MS/MS) study for the identification and characterization of in vivo metabolites of cisplatin in rat kidney cancer tissues: online hydrogen/deuterium (H/D) exchange study, *PLoS. One* 10 (8) (2015) e0134027.
- [15] N. Kumar, N. Goel, Phenolic acids: natural versatile molecules with promising therapeutic applications, *Biotechnol. Rep.* 24 (2019) e00370.
- [16] D. Parmar, E. Sugiono, S. Raja, M. Rueping, Complete field guide to asymmetric BINOL-phosphate derived Brønsted acid and metal catalysis: history and classification by mode of activation; Brønsted acidity, hydrogen bonding, ion pairing, and metal phosphates, *Chem. Rev.* 114 (18) (2014) 9047–9153.
- [17] U.H. Verkerk, P. Kébarle, Ion-ion and ion-molecule reactions at the surface of proteins produced by nanospray. Information on the number of acidic residues and control of the number of ionized acidic and basic residues, *J. Am. Soc. Mass Spectrom.* 16 (2005) 1325–1341.
- [18] W.-U. Palm, N. Schmidt, M. Stahn, S. Grimme, A kinetic study of the photolysis of sulfamethoxazole with special emphasis on the photoisomer, *Photochem. Photobiol. Sci.* 22 (3) (2023) 615–630.
- [19] S. Mezghich, F. Ayari, M. Sarakha, Photochemical degradation of the antimicrobial sulfamethoxazole upon solar light excitation: kinetics and elucidation of byproducts using LC/ESI+/MS2 technique, *Mass Spectrom. Purif. Tech.* 3 (2017) 1000118.
- [20] R.I.D. Birkler, N.B. Støttrup, S. Hermansson, T.T. Nielsen, N. Gregersen, H. E. Bøtker, M.F. Andreassen, M. Johannsen, A UPLC–MS/MS application for profiling of intermediary energy metabolites in microdialysis samples—A method for high-throughput, *J. Pharm. Biomed. Anal.* 53 (4) (2010) 983–990.
- [21] T. Cajka, J. Hricko, L. Rudl Kulhava, M. Paucova, M. Novakova, O. Kuda, Optimization of mobile phase modifiers for fast LC-MS-based untargeted metabolomics and lipidomics, *Int. J. Mol. Sci.* 24 (3) (2023) 1987.
- [22] T.H. Walter, B.A. Alden, J.A. Field, N.L. Lawrence, D.L. Osterman, A.V. Patel, M. A. DeLoffi, Characterization of a highly stable mixed-mode reversed-phase/weak anion-exchange stationary phase based on hybrid organic/inorganic particles, *J. Sep. Sci.* 44 (5) (2021) 1005–1014.
- [23] M.F. Medellín-Martínez, I. Luna-Zavala, M. Martínez-Delgado, J.T. Pérez-Urizar, J. A. Ramírez-Telles, O. Patiño-Rodríguez, Sensitive assay of clenbuterol residues in beef by ultra-high performance liquid chromatography coupled with mass spectrometry (UPLC-MS/MS) and solid-phase extraction, *Food Anal. Methods* 11 (2018) 2561–2568.
- [24] Z. Liu, Y. Shao, G. Zhu, X. Wang, Y. Chai, L. Wang, Gas phase reaction between chromones and solvent in an electrospray ionization source, *J. Mass Spectrom.* 54 (1) (2019) 66–72.
- [25] O. Núñez, P. Lucci, Applications and uses of formic acid in liquid chromatography–mass spectrometry analysis, *Adv. Chem. Res.* 20 (2013) 71–86.
- [26] H. Abdallah, C. Arnaudguilhem, R. Lobinski, F. Jaber, A multi-residue analysis of sulphonamides in edible animal tissues using QuEChERS extraction and HPLC-MS/MS, *Anal. Methods* 7 (4) (2015) 1549–1557.
- [27] A.H. Shendy, M.A. Al-Ghobashy, S.A.G. Alla, H.M. Lotfy, Development and validation of a modified QuEChERS protocol coupled to LC-MS/MS for simultaneous determination of multi-class antibiotic residues in honey, *Food Chem.* 190 (2016) 982–989.
- [28] N. Zhang, W. Deng, Y. Li, Y. Ma, Y. Liu, X. Li, H. Wang, Formic acid of ppm enhances LC-MS/MS detection of UV irradiation-induced DNA dimeric photoproducts, *Anal. Chem.* 92 (1) (2019) 1197–1204.
- [29] J. Chen, G.-G. Ying, W.-J. Deng, Antibiotic residues in food: extraction, analysis, and human health concerns, *J. Agric. Food Chem.* 67 (27) (2019) 7569–7586.
- [30] H.-N. Jung, D.-H. Park, K.-H. Yoo, H.-J. Cho, J.-H. Shim, H.-C. Shin, A. Abd El-Aty, Simultaneous quantification of 12 veterinary drug residues in fishery products using liquid chromatography–tandem mass spectrometry, *Food Chem.* 348 (2021) 129105.
- [31] J.A. Park, A. Abd El-Aty, W. Zheng, S.K. Kim, J.m. Choi, A. Hacımuftioğlu, J. H. Shim, H.C. Shin, Determination of metoserate, buquinolate, and diclofenac in pork, eggs, and milk using liquid chromatography–tandem mass spectrometry, *Biomed. Chromatogr.* 32 (6) (2018) e4215.
- [32] C. Mattos, D. Ringe, Proteins in organic solvents, *Curr. Opin. Struct. Biol.* 11 (6) (2001) 761–764.
- [33] A. Kaufmann, P. Butcher, K. Maden, S. Walker, M. Widmer, Multi-residue quantification of veterinary drugs in milk with a novel extraction and cleanup technique: salting out supported liquid extraction (SOSLE), *Anal. Chim. Acta* 820 (2014) 56–68.
- [34] A.M. Escorsim, G. da Rocha, J.V. Vargas, A.B. Mariano, L.P. Ramos, M.L. Corazza, C.S. Cordeiro, Extraction of Acutodesmus obliquus lipids using a mixture of ethanol and hexane as solvent, *Biomass Bioenergy* 108 (2018) 470–478.
- [35] R. Halim, R. Harun, M.K. Danquah, P.A. Webley, Microalgal cell disruption for biofuel development, *Appl. Energy* 91 (1) (2012) 116–121.
- [36] F.J. Baumgardt, A. Zandoná Filho, M.V. Brandalize, D.C. da Costa, N.R. Antoniosi Filho, P.C. Abreu, M.L. Corazza, L.P. Ramos, Lipid content and fatty acid profile of Nannochloropsis oculata before and after extraction with conventional solvents and/or compressed fluids, *J. Supercrit. Fluids* 108 (2016) 89–95.
- [37] S.-H. Chang, Y.-H. Lai, C.-N. Huang, G.-J. Peng, C.-D. Liao, Y.-M. Kao, S.-H. Tseng, D.-Y. Wang, Multi-residue analysis using liquid chromatography tandem mass spectrometry for detection of 20 coccidiostats in poultry, livestock, and aquatic tissues, *J. Food Drug Anal.* 27 (3) (2019) 703–716.
- [38] D. Zhang, Z.-Y. Park, J.-A. Park, S.-K. Kim, D. Jeong, S.-H. Cho, J.-H. Shim, J.-S. Kim, A. Abd El-Aty, H.-C. Shin, A combined liquid chromatography–triple-quadrupole mass spectrometry method for the residual detection of veterinary drugs in porcine muscle, milk, and eggs, *Environ. Monit. Assess.* 188 (2016) 1–16.
- [39] M. Schmid, F. Guihéneuf, D.B. Stengel, Evaluation of food grade solvents for lipid extraction and impact of storage temperature on fatty acid composition of edible seaweeds *Laminaria digitata* (Phaeophyceae) and *Palmaria palmata* (Rhodophyta), *Food Chem.* 208 (2016) 161–168.
- [40] Q. Guo, Y. Yang, J. Li, B. Shao, J. Zhang, Screening for plant toxins in honey and herbal beverage by ultrahigh-performance liquid chromatography-ion mobility-quadrupole time of flight mass spectrometry, *Am. J. Anal. Chem.* 13 (3) (2022) 108–134.
- [41] U. Koesukwiwat, S.J. Lehotay, K. Mastovska, K.J. Dorweiler, N. Leepipatpiboon, Extension of the QuEChERS method for pesticide residues in cereals to flaxseeds, peanuts, and doughs, *J. Agric. Food Chem.* 58 (10) (2010) 5950–5958.
- [42] X. Wu, T. Li, H. Feng, Y. Xie, F. Liu, K. Tong, C. Fan, Y. Liu, H. Chen, Multi-residue analysis of 206 pesticides in grass forage by the one-step quick, easy, cheap, effective, rugged, and safe method combined with ultrahigh-performance liquid chromatography quadrupole orbitrap mass spectrometry, *J. Sep. Sci.* 45 (14) (2022) 2520–2528.
- [43] F.J. Schenck, P. Callery, P.M. Gannett, J.R. Daft, S.J. Lehotay, Comparison of magnesium sulfate and sodium sulfate for removal of water from pesticide extracts of foods, *J. AOAC Int.* 85 (5) (2002) 1177–1180.
- [44] E. Chambers, D.M. Wagrowski-Diehl, Z. Lu, J.R. Mazzeo, Systematic and comprehensive strategy for reducing matrix effects in LC/MS/MS analyses, *J. Chromatogr. B* 852 (1–2) (2007) 22–34.
- [45] H. Stahnke, T. Reemtsma, L. Alder, Compensation of matrix effects by postcolumn infusion of a monitor substance in multiresidue analysis with LC–MS/MS, *Anal. Chem.* 81 (6) (2009) 2185–2192.
- [46] M. Urpi-Sardà, O. Jauregui, R.M. Lamuela-Raventos, W. Jaeger, M. Miksits, M.-I. Covas, C. Andres-Lacueva, Uptake of diet resveratrol into the human low-density lipoprotein. Identification and quantification of resveratrol metabolites by liquid chromatography coupled with tandem mass spectrometry, *Anal. Chem.* 77 (10) (2005) 3149–3155.
- [47] N. Van Hoof, D. Courtheyn, J.P. Antignac, M. Van de Wiele, S. Poelmans, H. Noppe, H. De Brabander, Multi-residue liquid chromatography/tandem mass spectrometric analysis of beta-agonists in urine using molecular imprinted polymers, *Rapid Commun. Mass Spectrometry* 19 (19) (2005) 2801–2808.
- [48] J. Bhatt, G. Subbaiah, S. Singh, Liquid chromatography/tandem mass spectrometry method for simultaneous determination of risperidone and its active metabolite 9-hydroxyrisperidone in human plasma, *Rapid Commun. Mass Spectrom.* 20 (14) (2006) 2109–2114.

- [49] L. Sabatini, A. Barbieri, M. Tosi, A. Roda, F.S. Violante, A method for routine quantitation of urinary 8-hydroxy-2'-deoxyguanosine based on solid-phase extraction and micro-high-performance liquid chromatography/electrospray ionization tandem mass spectrometry, *Rapid Commun. Mass Spectrometry* 19 (2) (2005) 147–152.
- [50] E. Stokvis, H. Rosing, J.H. Beijnen, Stable isotopically labeled internal standards in quantitative bioanalysis using liquid chromatography/mass spectrometry: necessity or not? *Rapid Commun. Mass Spectrometry* 19 (3) (2005) 401–407.
- [51] A. Kloepper, J.B. Quintana, T. Reemtsma, Operational options to reduce matrix effects in liquid chromatography–electrospray ionisation-mass spectrometry analysis of aqueous environmental samples, *J. Chromatogr. A* 1067 (1–2) (2005) 153–160.
- [52] M. Petrović, M.D. Hernando, M.S. Díaz-Cruz, D. Barceló, Liquid chromatography–tandem mass spectrometry for the analysis of pharmaceutical residues in environmental samples: a review, *J. Chromatogr. A* 1067 (1–2) (2005) 1–14.
- [53] J.B. Quintana, T. Reemtsma, Sensitive determination of acidic drugs and triclosan in surface and wastewater by ion-pair reverse-phase liquid chromatography/tandem mass spectrometry, *Rapid Commun. Mass Spectrometry* 18 (7) (2004) 765–774.
- [54] Z. Hu, P. Qi, N. Wang, Q.-Q. Zhou, Z.-H. Lin, Y.-Z. Chen, X.-W. Mao, J.-J. Jiang, C. Li, Simultaneous determination of multiclass illegal dyes with different acidic–basic properties in foodstuffs by LC-MS/MS via polarity switching mode, *Food Chem.* 309 (2020) 125745.
- [55] Codex Alimentarius Commission(2013). CCRVDF Codex Committee on Residues of Veterinary Drugs in Foods (2012) Guidelines for the design and implementation of national regulatory food safety assurance programme associated with the use of veterinary drugs in food producing animals (CAC/GL 71-2009).
- [56] E. Kim, H. Park, S. Park, J. Choi, H.J. Yoon, J.-H. Kim, Simultaneous determination of multi-class veterinary drugs in fishery products with liquid chromatography–tandem mass spectrometry, *Appl. Biol. Chem.* 64 (1) (2021) 40.
- [57] S.S. Rha, Y.J. Yang, W.H. Kim, Y.B. Jin, K.I. Park, H.-J. Lee, Establishment of the Maximum Residual Limit in the Milk of Dairy Cows Injected Intramuscularly with Prednisolone, *Vet. Sci.* 10 (10) (2023) 614.
- [58] Korean Ministry of Food and Drug Safety(2015). Applicability of Veterinary Drug MRLs for food in general. <https://residue.foodsafetykorea.go.kr/vd/mrl>. (Accessed 11 May 2024).
- [59] Government of Canada (2024). Maximum Residue Limits (MRLs) for Veterinary Drugs in Foods. <https://www.canada.ca/en/health-canada/services/drugs-health-products/veterinary-drugs/maximum-residue-limits-mrls/list-maximum-residue-limits-mrls-veterinary-drugs-foods.html>. (Accessed 13 May 2024).
- [60] US Food and Drug Administration(2023). CFR - Code of Federal Regulations Title 21, www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfCFR/CFRSearch.cfm?CFRPart=556&showFR=1. (Accessed 11 May 2024).
- [61] M. Kim, N. Kim, H. Kwon, S. Ha, H. Kim, J. Kim, Development of a simultaneous multi-residue analysis for screening and confirmation of 7 veterinary drugs in bovine milk by LC-MS/MS, *J. Prev. Vet. Med.* 43 (2019) 68–73.
- [62] H. DiMaio Knych, E. Steffey, J. Deuel, R. Shepard, S. Stanley, Pharmacokinetics of yohimbine following intravenous administration to horses, *J. Vet. Pharmacol. Ther.* 34 (1) (2011) 58–63.
- [63] D. Lucas, J. Neal-Kababick, J. Zweigenbaum, Characterization and quantitation of yohimbine and its analogs in botanicals and dietary supplements using LC/QTOF-MS and LC/QQQ-MS for determination of the presence of bark extract and yohimbine adulteration, *J. AOAC Int.* 98 (2) (2015) 330–335.
- [64] G.O. Quiroz-Delfi, C.V. Rider, S.S. Ferguson, A.K. Jarmusch, G.A. Mueller, Non-targeted chemical analysis of consumer botanical products labeled as blue cohosh (*Caulophyllum thalictroides*), goldenseal (*Hydrastis canadensis*), or yohimbe bark (*Pausinystalia yohimbe*) by NMR and MS, *Anal. Bioanal. Chem.* 416 (1) (2024) 175–189.
- [65] J. Sun, P. Chen, Chromatographic fingerprint analysis of yohimbe bark and related dietary supplements using UHPLC/UV/MS, *J. Pharm. Biomed. Anal.* 61 (2012) 142–149.
- [66] T. Nakajima, H. Hayashi, T. Sasamoto, M. Kanda, T. Kusano, Y. Matsushima, S. Kanai, K. Takeba, T. Nagayama, Determination of nequinat and buquinolone in livestock products using liquid chromatography–tandem mass spectrometry. *Shokuhin eiseigaku zasshi, J. Food Hyg. Soc. Japan* 52 (3) (2011) 178–182.
- [67] Y.S. Jung, D.-B. Kim, T.G. Nam, D. Seo, M. Yoo, Identification and quantification of multi-class veterinary drugs and their metabolites in beef using LC–MS/MS, *Food Chem.* 382 (2022) 132313.
- [68] J.A. Park, D. Zhang, D.S. Kim, S.K. Kim, K.S. Cho, D. Jeong, J.H. Shim, J.S. Kim, A. Abd El-Aty, H.C. Shin, Single-step multiresidue determination of ten multiclass veterinary drugs in pork, milk, and eggs using liquid chromatography with tandem mass spectrometry, *J. Sep. Sci.* 38 (16) (2015) 2772–2780.
- [69] A.P. Pais, R.S.B. Franco, A.F.B. da Silva, O. de Freitas, J.A.R. Paschoal, LC-MS/MS methods for sulfadimethoxine and ormetoprim analysis in feed and fish fillet and a leaching study for feed after alternative procedures for the incorporation of drugs, *Food Addit. Contam. Part A* 34 (4) (2017) 501–508.
- [70] H. Moon, A.-j Nam, K.E. Muambo, J.-E. Oh, Simultaneous multi-residue analytical method for anesthetics and sedatives in seafood samples by LC-ESI/MS/MS, *Food Chem.* 404 (2023) 134157.